

Product datasheet for **MR215347**

Zfp568 (NM_001167873) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Zfp568 (NM_001167873) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Zfp568
Synonyms:	C80731; chato; Gm1691; Znf568
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>MR215347 representing NM_001167873
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGAGCGCTTGTCCAGATGGCAGGAAGAAGAGCTTGGTGTGCCGAGGATTCTGTACCTCGCCAGGAAG
 AAGAAGACAGGACCAGGCCTTCAAAAACAGTGACATTCAGGATGTGGCTGTGGACCTCACACAGGAGGA
 ATGGCAGCAAATGAAGCCTGCTCAGAGGCTTTGTACCGAGATGTGATGTGGAGACCTACAGCAACCTA
 GTAACCGTGGGTGTCAAGTCACCAAACAGATGTGATCTTCAAGTTGGAACAAGCAGAGGAGCCATGGG
 TACTGGAGGAAGAGATGTTTGGAGACGCTCTCCAGCCGCAAGAGGAAGAATGAAGAGTTTTGCCTTTAA
 GGACATGGCCAAGGATCTGAGGTTTCAGGATGTGGTCATTTACTTCTCTCTGGAGGAGTGGGAGTGCTTG
 CGCCATTTCCACAGGAATCTGTACCGAGCTGTGATGTTGGATAACTACAGCAACCTGCTGTCACTGAGTC
 TTGCTGATACTAAGCCACGGGTGGTCTCCCTTCTAGAACAGGGCAAAGAGCCCTGGATGGTTATGAGGAA
 TGAGACAAAAATTTGGCATCCTGATTGGGTGTCTAGGACAGAAGCCAAGGACAGTTCCAAAATTAATAACA
 TTGCAAGAAAAAGATGGCAAAAAAACACACATGCCCCACTCTTGAGGACAGCAAGACCAGAGGGGACAGAG
 AAGTCACACGTGAGCTGGAGGGTCAACAGTCCATCAGGAGGGACACTTGAGGCAAGCTGCAGTAACCTC
 TGTCGAGAGGCCTGATTCTGTCCAGTGCACAGCCCATAGGGAGGCTCACCTGGAGGAAAACCTGCAGC
 TCTGAGAAAAGTCAGAAAACGTCCTGTGCCAGCCACCTCCCATTTGAACGTGAGCAACTTCACAGTAAGG
 CAAAGGCCAGTGAGCATGCACAGCATGGCAAGGTCTTTAATAGCTGCACGTGACAGACGGCTGTGCATCC
 GAGACCCAGGAAAGCAGAAAAGACAGTGAGAGGAAGAAAAGTGCCTTGCTGGTGGTCTGACACCTCG
 AAACCTCAGAGCGCTCAGGGGAGTGAGAGGCCACATAAGTGTAAAGGAATGCGGGAAGGCCTTCCACACTC
 CCTCCAGCTCAGTCACCATCAAAAGCTGCATGTGGGTGAGAAGCCCTACAAGTGCCAAGAGTGTGAAAA
 GGCTTTCCCGTCTAACGCCAGCTTAGTCTTACCATAGAGTTCACACTGATGAGAAGTGCTTTGAATGC
 AAGGAGTGCGGGAAGCCTTCATGCGCCCATCACACCTGCTCCGGCACCGCGCATCCACACTGGGGAGA
 AGCCTCACAAGTGTAAAGGAGTGCGGGAAGGCCTTCCGCTACGACACACAGCTCAGCCTGCACCTGCTCAC
 GCATGCTGGGGCACGGCGCTTTGAGTGTAAAGACTGTGACAAGGTGTATAGCTGTGCCTCCAGCTCGCT
 CTCCACCAGATGTCACACACTGGGGAGAAGCCCCACAAGTGTAAAGGAGTGTTGGAAAGGCTTCATCTCTG
 ACTCACACCTTCTCGCCATCAGAGTGTTCACACTGGCGAGACGCCCTATAAGTGTAAAGGAGTGTGGCAA
 AGGCTTCCGTCTGGTTCAGAACTTGCCCGGCATCAGAGAGCCATTCTGGCGATAAACCTACAAGTGT
 AAGGAATGTGAAAGTCCTTCACTTGACACACAGAGCTCTTTAGACATCAGAAAGTTCACACTGGTGATA
 GACCCCAAAATGTAAAGGAGTGTGGGAAGGCCTTCATCCGAAGATCAGAGCTGACACATCATGAAAGAAG
 CCACAGCGGTGAGAAACCTATGAGTGCAAAGAGTGTGGAAGACCTTCGGCCGAGGCTCAGAGCTCAGT
 CGACACCAGAAAAATCCACACCGGTGAGAAGCCCTATAAGTGCCAGCAATGTGGGAAGGCCTTTATTCGTG
 GCTCTACCTTACCAACATCAGAGGATTCATACAGGACGGAGGAGTGAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >MR215347 representing NM_001167873
Red=Cloning site Green=Tags(s)

MERLSQMAGRAWCAEDSVPRQEEEDRTRPSKTVTFKDVAVDLTQEEWQQMKPAQRALYRDVMLETYSNL
 VTVCQVTKPDVIFKLEQAEPPWVLEEMFWRRSPAARGRMKSFAFKDMAKDLRFEDVVIYFSLEEWEC
 RHSHRNLYRAVMDNYSNLLSLSLADTKPRVVSLLLEQGKEPMMVMRNETKIWHPDWVSRTAKDSSKI
 LQEKMAKKHTCPTLEDSKTRGDREVTRELEGQQVHQEGHLRQAAVTSVERPDSVQCTAHREAHPGGKPCS
 SEKSQKTSKLCQPPPIEREQLHSAKASEHAQHGVFNSTSDTAVHPRPQESRKDSERKKSALAGGPDTS
 KPQSAQGSERPHKCECGKAFHTPSQLSHHQKLHVGEKPYKCQCEGKAFPSNAQLSLHHRVHTDEKCFEC
 KECGKAFMRPSHLLRHQRIHTGEKPHKCECGKAFRYDTQLSLHLLTHAGARRFECKDCDKVYSCASQLA
 LHQMSHTGEKPHKCECGKGFISDHLRLHQSHTGETPYKCECGKGFRRGSELARHQAHSQDKPYKC
 KECGKSFTCTTELFRHQVHTGDRPHKCECGKAFIRSELTHHERSHSGEKPYECKEKGTFGRGSELS
 RHQKIHTGEKPYKCQCGKAFIRGSHLTQHRIHTGRSE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Cloning Scheme:



ACCN: NM_001167873

ORF Size: 2010 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001167873.1](#), [NP_001161345.1](#)

RefSeq Size: 7199 bp

RefSeq ORF: 2013 bp

Locus ID: 243905


UniProt ID: [E9PY11](#)

Cytogenetics: 7 B1

MW: 77.4 kDa

Gene Summary: Has transcriptional repression activity, partially through the recruitment of the corepressor TRIM28 but has also repression activity independently of this interaction (PubMed:22110054, PubMed:27658112). Essential during embryonic development, where it acts as direct repressor of IGF2-P0, placental-specific transcript of IGF2, in early development and regulates convergent extension movements required for axis elongation and tissue morphogenesis in all germ layers (PubMed:18701545, PubMed:22110054, PubMed:28522536). Also important for normal morphogenesis of extraembryonic tissues including the yolk sac, extraembryonic mesoderm and placenta (PubMed:18701545, PubMed:21094155). May enhance proliferation or maintenance of neural stem cells (PubMed:23071813).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR215347

