

Product datasheet for **MR215358L4V**

Scel (NM_022886) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Scel (NM_022886) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Scel
Synonyms:	9230114I02Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_022886
ORF Size:	1959 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215358).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_022886.2 , NP_075024.2
RefSeq Size:	3053 bp
RefSeq ORF:	1959 bp
Locus ID:	64929
UniProt ID:	Q9EQG3
Cytogenetics:	14 E2.3
Gene Summary:	May function in the assembly or regulation of proteins in the cornified envelope. The LIM domain may be involved in homotypic or heterotypic associations and may function to localize sciellin to the cornified envelope (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)