

Product datasheet for **MR215365A1V**

Mouse Ccl26 (NM_001013412) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Ccl26 (NM_001013412) AAV Particle
Tag: Myc-DDK
Symbol: Ccl26
Synonyms: Ccl26l
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR215365 representing NM_001013412
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTTCTTCGATTTGGGTCTCCTTGTTCCTGGCCATCTTCCTCAGTGTCCAGCTTGGTGTGCCACAT
GTGGGAGCAGCATCGCTATGCTCTGCTGCCCTAATTCAGCTACTATGTGATCCCATGGAGCTGGGTGTA
CAGCTATAAGTTCACCGACAAGAGCTGCACCACTGACGGTGTGATATTCTTTACAAAAACAGGTAAGCAA
TTCTGTGTCCAGCCAGGGGCCAAATGGGTGCAGAGATTCATCTCTCTGGTGAACACCAGGAACCAATTTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR215365 representing NM_001013412
Red=Cloning site Green=Tags(s)

MFFDLGLLVLLAIFLSVQLGVATCGSSIAMSCCPNFSYVYVIPWSWVYSYKFTDKSCTSDGVIFFTKTGKQ
FCVQPGAKWVQRFISLVNTRNHL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_001013412
ORF Size: 279 bp
Buffer: PBS with 0.001% Pluronic F68



[View online »](#)

| | |
|----------------------|---|
| Stability: | AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus. |
| RefSeq: | <u>NM_001013412.2</u> , <u>NP_001013430.2</u> |
| RefSeq Size: | 282 bp |
| RefSeq ORF: | 282 bp |
| Locus ID: | 541307 |
| Cytogenetics: | 5 G2 |
| MW: | 10.5 kDa |