

## Product datasheet for **MR215380L3V**

### Scarf1 (NM\_001004157) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Scarf1 (NM_001004157) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Scarf1
Synonyms:	AA986099; mKIAA0149; SREC; SREC-I
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001004157
ORF Size:	2460 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215380).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001004157.2</a> , <a href="#">NP_001004157.2</a>
RefSeq Size:	2884 bp
RefSeq ORF:	2463 bp
Locus ID:	380713
UniProt ID:	<a href="#">Q5ND28</a>
Cytogenetics:	11 B5
Gene Summary:	Mediates the binding and degradation of acetylated low density lipoprotein (Ac-LDL). Mediates heterophilic interactions, suggesting a function as adhesion protein (By similarity). Plays a role in the regulation of neurite-like outgrowth.[UniProtKB/Swiss-Prot Function]



[View online »](#)