

## Product datasheet for **MR215391L3V**

### Scnm1 (NM\_001163573) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Scnm1 (NM_001163573) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Scnm1
Synonyms:	3110001117Rik; Scnm1-ps
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001163573
ORF Size:	492 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215391).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001163573.1</a> , <a href="#">NP_001157045.1</a>
RefSeq Size:	693 bp
RefSeq ORF:	495 bp
Locus ID:	69269
UniProt ID:	<a href="#">Q8K136</a>
Cytogenetics:	3 40.74 cM



[View online »](#)

**Gene Summary:**

Mutations in the voltage-gated sodium channel gene Scn8a lead to neurological problems in mice. For one particular mutation, Scn8amedJ, mice live to adulthood but have tremors and muscle weakness, among other problems, in all strains except those derived from C57BL6 mice. In these strains, the product of the Scnm1 gene (229 aa) partially overcomes the effects of the Scn8amedJ mutation. However, in C57BL6-derived mice, a one nt change in the penultimate exon creates a premature stop codon, truncating the Scnm1 protein at 186 aa. This truncated protein lacks the ability to overcome the effects of the Scn8amedJ mutation, and these mice suffer paralysis and juvenile death. [provided by RefSeq, Jul 2009]