

## Product datasheet for **MR215393L3V**

### **Pcsk6 (NM\_011048) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Pcsk6 (NM_011048) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pcsk6
Synonyms:	b2b2830Clo; C86343; Pace4; SPC4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_011048
ORF Size:	2877 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215393).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_011048.1</a> , <a href="#">NP_035178.1</a>
RefSeq Size:	4298 bp
RefSeq ORF:	2880 bp
Locus ID:	18553
Cytogenetics:	7 35.36 cM



[View online »](#)