

## Product datasheet for **MR215397L3V**

### Ctsw (NM\_009985) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ctsw (NM_009985) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ctsw
Synonyms:	lym; lymphopain
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_009985
ORF Size:	1113 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215397).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_009985.5</a> , <a href="#">NP_034115.2</a>
RefSeq Size:	1438 bp
RefSeq ORF:	1116 bp
Locus ID:	13041
UniProt ID:	<a href="#">P56203</a>
Cytogenetics:	19 A



[View online »](#)

**Gene Summary:**

This gene encodes a member of the peptidase C1 (papain) family of cysteine proteases. The encoded preproprotein is proteolytically processed to generate a mature protein product. Expression of the encoded protein is upregulated following lymphocyte activation. Data from a human cell line suggests that the encoded enzyme may be important for viral entry into host cells. [provided by RefSeq, Aug 2015]