

## Product datasheet for **MR215424**

### **Tmem176b (NM\_001164208) Mouse Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Tmem176b (NM\_001164208) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Tmem176b  
**Synonyms:** 1810009M01Rik; Clast1; Lr8  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >MR215424 representing NM\_001164208  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGTCCAGAGCACAGTGACTGTGAATGGAGTCAAAGTGGCCTCCACACATCCCCAGTCCGCTCACATCA  
GCATCCACATCCACAAAAATCTGCTTTGGAGCAGCTGCTGGGAGCTGTGGGCTCCCTGAAGAAGTTTCT  
CTCCTGGCCTCAGGCCAGAATCCACTATGGGCAGCTGTCTCTAGGGGTGACCCAGATATTGCTGGGGCTT  
GTGAGCTGTGCTCTTGGAGTATGTCTCTACTTCGGGCCTTGACTGAGCTGTGCGCCTTTGGCTGTGCTT  
TCTGGTCGGGTCTGTGGCAATTCTAGCTGGAGTTGGTACTATTGTCCATGAGAAGAGGCAGGGAAAAC  
GTCCGGCCAGGTATCATGTCTCCTCCTGGCTTGCAATGCCACAGCTGCAGCTGTACGTTCTGGGT  
GTGAACAGCTTAATCAGGCAAACAGTGTTCCTACTACGTCGAGATCTTTCCACATGCAACCCCTTGC  
AATCAAGCATGGACCCTGGGTATGGGACAGTGCATACAGTGTGATTCAGACTGGAAGACAGAGAGGTG  
CAGAGAGTACCTGAACATGATGATGAACCTGTTCCCTAGCGTTCTGCATCATGCTCACAGTGGTCTGTATC  
CTGGAGATTGTTGTCTGTGGCTTCCCTGGGACTAAGTCTCCGAAGTATGTATGGCCGGAGCTCCCAGG  
CCCTGAATGAAGAAGAGTCAGAGAGGAAGCTTCTTGATGGCATCCAGCACCTGCGTCCCCGGCCAAAGA  
GAAGATCCCTGCTATCCTG

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >MR215424 representing NM\_001164208  
 Red=Cloning site Green=Tags(s)

MVQSTVTVNGVKVASTHPQSAHISIHQKSALEQLLGAVGSLKKFLSWPQARIHYGQLSLGVTQILLGL  
 VSCALGVCLYFGPWTELCAFGCAFWSGSVAILAGVGTIVHEKRQKLSGQVSCLLLLACIATAAAATVLG  
 VNLIHQTSVPYYVEIFSTCNPLQSSMDPGYGTVRYSDSDWKTERCREYLNMMNMLFLAFCIMLTVVCI  
 LEIVVSVASLGLSLRSMYGRSSQALNEEESERKLLDGHAPAPASPAKEKIPAIL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF

**ACCN:** NM\_001164208

**ORF Size:** 789 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001164208.1](#), [NP\\_001157680.1](#)

**RefSeq Size:** 1069 bp

**RefSeq ORF:** 792 bp

**Locus ID:** 65963


**UniProt ID:** [Q9R1Q6](#)

**Cytogenetics:** 6 23.75 cM

**MW:** 28.4 kDa

**Gene Summary:** May play a role in the process of maturation of dendritic cells (By similarity). Required for the development of cerebellar granule cells.[UniProtKB/Swiss-Prot Function]

### Product images:



Circular map for MR215424