

Product datasheet for **MR215431**

Zfp568 (NM_001033355) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	Zfp568
Synonyms:	C80731; chato; Gm1691; Znf568
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



ORF Nucleotide Sequence: >MR215431 representing NM_001033355
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGAGCGCTTGTCAGATGGCAGGAAGAAGAGCTTGGTGTGCCGAGGATTCTGTACCTCGCCAGGAAG
 AAGAAGACAGGACCAGGCTTCAAAAACAGTGACATTCAAGGATGTGGCTGTGGACCTCACACAGGAGGA
 ATGGCAGCAAAATGAAGCCTGCTCAGAGGGCTTTGTACCGAGATGTGATGCTGGAGACCTACAGCAACCTA
 GTAACCGTGGGTTGTCAAGTCAACAAACAGATGTGATCTTCAAGTTGGAACAAGCAGAGGAGCCATGGG
 TACTGGAGGAAGAGATGTTTTGGAGACGCTCTCCAGAAGCCGCAAGAGGAAGAATGAAGAGTTTTGCCTT
 TAAGGACATGGCCAAGGATCTGAGGTTGAGGATGTGGTCATTTACTTCTCTGAGGAGTGGGAGTGC
 TTGCGCCATTCCCACAGGAATCTGTACCGAGCTGTGATGTTGGATAACTACAGCAACCTGCTGTCACTGA
 GTCTTGCTGATACTAAGCCACGGGTGGTCTCCCTTCTAGAACAGGGCAAGAGCCCTGGATGGTTATGAG
 GAATGAGACAAAAATTTGGCATCTGATTGGGTGTCTAGGACAGAAGCCAAGGACAGTTCAAAAATTA
 ACATTGCAAGAAAAGATGGCAAAAAACACACATGCCCACTCTTGAGGACAGCAAGACCAGAGGGGACA
 GAGAAGTCACACGTGAGCTGGAGGGTCAACAGGTCCATCAGGAGGACACTTGAGGCAAGCTGCAGTAAC
 CTCTGTGAGAGGCTGATTCTGTCCAGTGCACAGCCCATAGGAGGGCTCACCTGGAGGAAAACCTGC
 AGCTCTGAGAAAAGTCAGAAAACGTCCCTGTGCCAGCCACCTCCATTGAACGTGAGCAACTTCACAGTA
 AGGCAAGGGCCAGTGAGCATGCACAGCATGGCAAGGCTTTAATAGCTGCACGTGAGACAGGCTGTGCA
 TCCGAGACCCAGGAAAGCAGAAAAGACAGTGAGAGGAAGAAAAGTGCCCTTGCTGGTGGTCTGACACC
 TCGAAACCTCAGAGCGCTCAGGGGAGTGAGAGGCCACATAAGTGTAAAGGATGCGGGAAGGCCTTCACA
 CTCCCTCCAGCTCAGTCACCATCAAAGCTGCATGTGGGTGAGAAGCCCTACAAGTGCCAAGAGTGTGG
 AAAGGCTTTCCCGTCTAACGCCAGCTTAGTCTTACCATAGAGTTTCACTGATGAGAAGTGCTTTGAA
 TGCAAGGAGTGCGGGAAGCCTTCATGCGCCCATCACACCTGCTCCGGCACCAGCGCATCCACACTGGGG
 AGAAGCCTCACAAGTGTAAAGAGTGCGGGAAGGCCTTCCGCTACGACACACAGCTCAGCTGCACCTGCT
 CACGCATGCTGGGACAGGCGCTTTGAGTGTAAAGACTGTGACAAGGTGTATAGCTGTGCTCCAGCTC
 GCTCTCCACAGATGTACACACTGGGGAGAAGCCCAAGTGTAAAGAGTGTTGGGAAAGGCTTCATCT
 CTGACTCACACCTTCTTCGCCATCAGAGTGTTCACACTGGCGAGACGCCCTATAAGTGTAAAGAGTGTGG
 CAAAGGCTTCCGTCGTGGTTCAGAACTTGGCCGGCATCAGAGAGCCATTCTGGCGATAAACCTACAAG
 TGTAAGGAATGTGAAAGTCTTCACTTGACCACAGAGCTCTTAGACATCAGAAAGTTCACACTGGTG
 ATAGACCCACAAATGTAAAGAGTGTGGGAAGGCCTTCATCCGAAGATCAGAGCTGACACATCATGAAAG
 AAGCCACAGCGGTGAGAAACCTATGAGTGCAAAGAGTGTGGAAAGACCTTCGGCCGAGGCTCAGAGCTC
 AGTCGACACCAGAAAATCCACACCGGTGAGAAGCCCTATAAGTGCCAGCAATGTGGGAAGGCCTTTATTC
 GTGGCTCTACCTTACCAACATCAGAGGATTCATACAGGACGGAGGAGTGAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >MR215431 representing NM_001033355
 Red=Cloning site Green=Tags(s)

MERLSQMAGRAWCAEDSVPRQEEEDRTRPSKTVTFKDVAVDLTQEEWQMKPAQRALYRDVMLETYSNL
 VTVGCQVTKPDVIFKLEQAEEPWVLEEEMFWRRSPEAARGRMKSF AFKDMAKDLRFEDVVIYFSLEEWEC
 LRHSHRNLYRAVMLDNYSNLLSLSLADTKPRVVSLLQKGEPWMVMRNETKIWHPDWVSRTAKDSSKIK
 TLQEKMAKKHTCPTLEDSKTRGDREVTRELEGQVHQEGHLRQAAVTSVERPDSVQCTAHREAHPPGKPC
 SSEKSQKTSKCPPIEREQLHSAKASEHAHQGVNSCTSDTAVHPRPQESRKDSERKKSALAGGPD
 SKPQSAQGSERPHKCECGKAFHTPSQLSHHQKLHVGEKPYKCECGKAFPSNAQLSLHHRVHTDEKCFE
 CKECGKAFMRP SHLLRHQRIHTGEKPHKCECGKAFRYDTQLSLHLLTHAGARRFECKDCDKVYSCASQL
 ALHQMSHTGEKPHKCECGKGFISD SHLLRHQSVHTGETPYKCECGKGFRRGSELARHQRAHSGDKPYK
 CKECGKSFTCTTELFRHQVHTGDRPHKCECGKAFIRRESELTHHERSHSGEKPYECKEKGKTFGRGSEL
 SRHQKIHTGEKPYKQCGKAFIRGSHLTQHRIHTGRRSE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001033355

ORF Size: 2013 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.


RefSeq: [NM_001033355.3](#), [NP_001028527.2](#)

RefSeq Size: 7195 bp

RefSeq ORF: 2016 bp
Locus ID: 243905
UniProt ID: [E9PYI1](#)
Cytogenetics: 7 B1
MW: 77.1 kDa

Gene Summary: Has transcriptional repression activity, partially through the recruitment of the corepressor TRIM28 but has also repression activity independently of this interaction (PubMed:22110054, PubMed:27658112). Essential during embryonic development, where it acts as direct repressor of IGF2-P0, placental-specific transcript of IGF2, in early development and regulates convergent extension movements required for axis elongation and tissue morphogenesis in all germ layers (PubMed:18701545, PubMed:22110054, PubMed:28522536). Also important for normal morphogenesis of extraembryonic tissues including the yolk sac, extraembryonic mesoderm and placenta (PubMed:18701545, PubMed:21094155). May enhance proliferation or maintenance of neural stem cells (PubMed:23071813).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR215431