

Product datasheet for **MR215471**

3632451O06Rik (NM_026142) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	3632451O06Rik (NM_026142) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	3632451O06Rik
Synonyms:	AU067705
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MR215471 representing NM_026142
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCTCCAGGACTCAATTACTGGAATTGTCAACTCCTTCAACCTGTTTTCCCATCCACTATGAGCAGAC
CAACTCTGATGCCACCTGTGTGGCTTTCTGTAGCATTCTTCCCTCACCTTGGCTACAGGATGTCAAGC
TTTCCAAAGGTGAAAAGGAGGGAGACAGCACAAGAATATGCAGAGAAAGAGCAGTCCCAGAAGATGAAC
ACTGATGACCAAGAAAACATCTCCTTTGCTCCAAAGTACATGCTGCAGCAGATGTCCTCTGAAGCTCCA
TGGTACTATCTGAAGGACCTTCAGAGATACCATTAATTAAGTATTCTCAGTCAACAAAGAGAGCCATCT
TCCAGGAGCTGGGCTTCTGCATCTACCAGCCCTGGTGTACTCTCCAGTGAAGCTGTGGTCTCAGCC
AGTGAACAAGAGCCTGGCCCCAGCCTGCTTGAAGAATGTCTTCTGAGCACAGCCTTTCAAAGTCAATGT
TAACTGTTGCTGTATCTCACCTGCTTCTGTAATCCTGACCAGGAGGGACCCTATAACAGTTTGAGTAC
CCAGCCTATTGTAGCAGCAGTACAGACGTAACACATGGTTCTCTGGATTATTTAGATAACCAAGTATTT
GCAGCTAAAAGTCAGGAAGCAGTTAGCTTAGGAAATTCACCTTCTTCTCTATAAAATACCAAAGAACCTG
AGATAATCAAAGCAGATGCAGCCATGGGAACACTGTCTGTTTCTGGTGTGATTCTACAGGAGACATGGA
GCCTGATAGGGAGAGACCCTCAGAGATGGCAGCTGATGATGGCCAGAGCACTACCACCAAGTATTTGGTG
ACAATCCCAATAATTTCTGACTACTGAGCCAACAGCTGGCAGTATTTTGGGGGATGCTAAGGTCACAG
TGAGTGTGAGCAGCAGCTGGCCAGTCTCTTCTATCTTCAATGAGGAGTGGGATGACACCAAAATTTGAGAG
CATAAGTCGGGGGAGGCCCCAGAGCCAGGAGACAATGCAGAGACTCAAATGAGAACAAGCCACCACAT
GGCACCTATGAGAGCTTTGAGGGGACTGAGGAGAGCCCTCTTCCACAGCAGTTTTAAAGGTGGCTCCAG
GGCACCTGGGAGGAGAACCAGCCTTGGGAACAGCCCTGGTGAATGCACTGGGGGATGAGAGGTCTCCAGT
CTTGACTCATCAGATTTCTTCACTCCCATGAGTCTGGCAGAAGACCCTGAAGTTTCTACCATGAAGCTG
TTTCAAGTGCAGGAGGCTTCAGAGCCTCCACCCAGGGGACAGGACCCAGCTCTCCTCAGAGACTGCAT
TTTCTACCTCTCAATATGAGTCTGTGCCACAGCAGGAAGCAGGAAATGTGTTAAAAGACATCACTCAAGA
GAGAAAAATGGCTACTCAAGCAATGAACACTACATCACCTGTGGTACACAAGAGCATATGGCCACCATC
GAGGTTCCAGAGGCAGTGGGGAGCCTGAAGAAGGCATGCCATCCCTGTCTCCTGTGCCTGCTGAAGTGG
CTGATGCTGAGCTGTCCAGAAGAGGGGAGTCTCTGGCCACTCCAGCCTCCACCACTGTTGTCCCTTTGTC
TCTCAAACCTTACCTTTCCATGGAAGATCTCATGGACACAATCACTGGACCAAGTGAAGAATTCATACCA
GTATTGGGCTCTCAATGGCACCCCTGCAATGACTGTGGAGGCTCCACAATTTCTTCAGCACTCCCTT
CTGAGGGAAGAAGTCCATCTATCAGCCGCCAAATACAGCTGCCTCCTATGGCTAGAACAACCTTGA
ATCTGAAGAGGTAGAAGATGATGAGGATGAAGAGGATGAAGAAGATGAAGAGGAAGAAGAAGATGAG
GAAGATGAAGAAGATGAGGAAGATAAAGAGACAGACTCCTTGATAAAGACTTCGATGGTGATGACTGAAC
CACCAGGTTTACTCTGCCTGGCATCACGTCCCAGGAACCTGACATACGTTTCAGGAAGCATGGACCTACT
AGAGGTAGCCACCTACCAGGTGCCAGAGACCATCGAGTGGGAACAGCAGAATCAAGGTCTGGTGAGAAGC
TGGATGGAAAAGCTAAAAGACAAGGCAGGCTACATGTCTGGAATGCTGGTGCCTGTTGGGGTAGGGATAG
CTGGAGCCTTGTTCATCTTGGGAGCGCTATACAGCATTAAAGTTATGAATCGCCGAAGAAGAAATGGATT
CAAAAGGCATAAAAGGAAGCAGAGAGAATTCAACAGTATGCAAGACAGAGTGATGCTGTTAGCTGATAGC
TCTGAAGACGAATTC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR215471 representing NM_026142
Red=Cloning site Green=Tags(s)


MLQDSITGIVNSFNLFFPSTMSRPTLMPTCVAFCSILFLTLATGCQAFPKVERRETAQEYAEKEQSQKMN
TDDQENISFAPKYMLQQMSSEAPMVLSEGPSEIPLIKVFSVNKESHLPGAGLLHPTSPGVYSSSEPVVSA
SEQEPGPSLLERMSSEHSLSKVMLTVAVSSPASLNPDQEGPYNSLSTQPIVAAVTDVTHGSLDYLDNQLF
AAKSQEAVSLGNSPSSSINTKEPEI IKADAAMGTTVVPGVDSTGMEPDRERPSEMAADDGQSTTTKYL
TIPNNFLTTEPTAGSILGDAKVTVSSTAGPVSSIFNEEWDDTKFESISRGRPPEPGDNAETQMRTKPPH
GTYESFEGTEESPSSTAVLKVAPGHLGGEPALGTALVTALGDERSPVLTHQISFTPMSLAEDPEVSTMKL
FPSAGGFRASTQGDRTQLSSETAFSTSQYESVPQQEAGNVLKDITQERKMATQAMNTTSPVVTQEHMATI
EVPRGSGEPEEGMPSPVPAEVADAELSRGESLATPASTTVVPLSLKLTSSMEDLMDTITGPSEEFIP
VLGSPMAPPAMTVEAPTISSALPSEGRTSPSISRPNNTAASYGLEQLESEEVEDDEDEDEDEDEDEDEDE
EDEDEEDKETSLSYKDFDGDTEPPGFTLPGITSQEPDIRSGSMDLLEVATYQVPETIEWEQNQGLVRS
WMEKLDKAGYMSGMLVPVGVGIAGALFILGALYSIKVMNRRRRNGFKRHKRQREFNSMQDRVMLLADS
SEDEF

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9004_b04.zip

Restriction Sites: Sgfl-Mlul

Cloning Scheme:



ACCN: NM_026142

ORF Size: 2325 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_026142.5](#)

RefSeq Size: 2748 bp

RefSeq ORF: 2328 bp


Locus ID: 67419

UniProt ID: [Q8BT18](#)

Cytogenetics: 14 C1

MW: 84.2 kDa

Product images:



Circular map for MR215471