

Product datasheet for **MR215471L1V**

3632451O06Rik (NM_026142) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	3632451O06Rik (NM_026142) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	3632451O06Rik
Synonyms:	AU067705
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_026142
ORF Size:	2325 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215471).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_026142.4
RefSeq Size:	2748 bp
RefSeq ORF:	2328 bp
Locus ID:	67419
UniProt ID:	Q8BT18
Cytogenetics:	14 C1



[View online »](#)