

Product datasheet for **MR215473L4V**

Adgrb1 (NM_174991) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Adgrb1 (NM_174991) Mouse Tagged ORF Clone Lentiviral Particle |
| Symbol: | Adgrb1 |
| Synonyms: | B830018M07Rik; Bai1; mKIAA4089; R75078 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_174991 |
| ORF Size: | 4746 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(MR215473). |
| OTI Disclaimer: | <p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p> |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_174991.3 , NP_778156.2 |
| RefSeq Size: | 6227 bp |
| RefSeq ORF: | 4749 bp |



[View online >](#)

Locus ID: 107831

UniProt ID: [Q3UHD1](#)

Cytogenetics: 15 D3

Gene Summary: Phosphatidylserine receptor which enhances the engulfment of apoptotic cells (PubMed:17960134). Also mediates the binding and engulfment of Gram-negative bacteria (PubMed:21245295, PubMed:26838550, PubMed:26838550). Stimulates production of reactive oxygen species by macrophages in response to Gram-negative bacteria, resulting in enhanced microbicidal macrophage activity (By similarity). In the gastric mucosa, required for recognition and engulfment of apoptotic gastric epithelial cells (By similarity). Promotes myoblast fusion (PubMed:23615608). Activates the Rho pathway in a G-protein-dependent manner (By similarity). Inhibits MDM2-mediated ubiquitination and degradation of DLG4/PSD95, promoting DLG4 stability and regulating synaptic plasticity (PubMed:25751059). Required for the formation of dendritic spines by ensuring the correct localization of PARD3 and TIAM1 (By similarity). Potent inhibitor of angiogenesis in brain and may play a significant role as a mediator of the p53/TP53 signal in suppression of glioblastoma (By similarity). [UniProtKB/Swiss-Prot Function]