

Product datasheet for **MR215474A1V**

Mouse Patell (NM_001128510) AAV Particle

Product data:

Product Type: AAV Particles

Tag: Myc-DDK

Symbol: Patell

Synonyms: 671003; Pate-P; Patell

Mammalian Cell None

Selection:

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >MR215474 representing NM_001128510
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGCATCGCC

ATGGGAAAGCACGCTTTCAGCTCCTGCTGGTCTGTCTTTGCTGGTGATGCCCTGCAAGCTTTGACAT
GTATCATGTGTAGGATGAATTCTCGGAGAATTTGTAGGGGAAAGAAGGTTGCTGTCAGGCTAGACC
TGGAGAGAAGTGTGCCTCGTTCATAACCTTAAAGATGGCAAAATTCAGTTTGAAGCCAGAGATGTGCT
GACCTTTGCTTCACTGGGACTGTTATGATTGGAGATAAGACAGTAAAAATGAACTGTTGCAACAACAAAT
CTTTCTGCAATAAACTA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR215474 representing NM_001128510
Red=Cloning site Green=Tags(s)

MGKHVLQLLLVLVSLVMSLQALTCITCDRMNSRRICEGKEGCCQARPGEKCSFITLKDGIQFGSQRCA
DLCFTGTVMIGDKTVKMNCNNKSFCKNL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse

Serotype: AAV-2

ACCN: NM_001128510

ORF Size: 297 bp



Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_001128510.1
RefSeq Size:	300 bp
RefSeq ORF:	300 bp
Locus ID:	671003
Cytogenetics:	9
MW:	10.9 kDa