

Product datasheet for MR215484L4

Nox3 (NM_198958) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Nox3 (NM_198958) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Nox3
Synonyms:	GP91-3; het; nmf25; nmf250
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215484).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_198958
ORF Size:	1704 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_198958.2 , NP_945196.2
RefSeq Size:	1792 bp
RefSeq ORF:	1707 bp
Locus ID:	224480
UniProt ID:	Q672J9
Cytogenetics:	17 2.05 cM
Gene Summary:	This gene encodes a member of the NOX family of NADPH oxidases. These enzymes catalyze the transfer of electrons from NADPH to molecular oxygen to produce superoxide and other reactive oxygen species (ROS). The ROS generated by family members have been implicated in numerous biological functions including host defense, posttranslational processing of proteins, cellular signaling, regulation of gene expression, and cell differentiation. The protein encoded by this gene is expressed predominantly in the inner ear and is involved in the biogenesis of otoconia, which are crystalline structures of the inner ear involved in the perception of gravity and linear acceleration. In mouse mutations of this gene lead to the absence of otoconia and vestibular dysfunction. [provided by RefSeq, Jun 2013]

Product images:



Circular map for MR215484L4