

Product datasheet for **MR215487**

Nolc1 (NM_053086) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Nolc1 (NM_053086) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Nolc1
Synonyms:	NOPP130; NOPP140; P130
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR215487 representing NM_053086
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCGGATACCGCTTGCGCCGCTGGTTCCAGCGACCTTTATCCCCTTGCTCAGATTTCTGCGGG
 ATAGCCAATCTCGGAGGTGGCCAGTAAATTTGCAAAAGCGACCGCGCTACACAGCAGGACGCCAATGC
 CTCGTCCTCTTGACATCTATAGCTTCTGGCTCAACAGGTCCACCAAAGCCCAAAGGTGAAGTTACAG
 TCAAATGGACAGTGACCAAGAAGGCTAAGAAAGAGACTTCATCCAGTGACAGCAGTGAGGACAGCAGTG
 AGGACGAGGACAAAAAGCCAGGGACTTCCACACAGAAGGCTGCCGCACAGGTCAAGCGAGCCAGTGT
 GCCTCAGCATGCTGAAAAGGCAGCAGCCAAAGCTTCCAGAGCAGCAGTGTGAAGAATCCAGTGAGGAA
 GAGGAAGAGGACAAAAAGAAAAGCCTGTCCAGAAGGCAGCTAAGCCCCAAGCCAAGGCAGTCAGACCTC
 CTGCGAAGAAGGCAGAGAGCTCTGAGTCGGACTCAGACTCGGATTCGGACTCCAGCTCAGAGGAAGAAAC
 ACCACAGACCCAGAAGCCAAAGGCAGCTGTGGCAGCAAAAGCTCAGACTAAAGCCGAAGCCAAACCAGGT
 ACACCAGCGAAAAGCACAGCCTAAGGTAGCCAATGGCAAAGCAGCCGCCAGCAGCAGCAGCAGCAGCAGCA
 GCGATGACTCAGAGGAAGAGAAGAGGCAGCTGCACCTCCCAAGAAGACTGTACCAAAAAAGCAAGTCTGT
 GGCCAAGGCCAGTGAAGTAGCTGCCGCCCCACCCAGAAGAGCTCCAGCAGTGAGGATCTTCCAGT
 GAAGAGGAGGAGGACAGAGACAACCCATGAAGAAAAAGCAGGTCCCTACAGTTCAGTTCACCCACCT
 CTGTTCTTTACCAAAGAAGTCCCCGGGAACCCAGGCTCCAAAGAAAGCTGCTGCGCAGACACAGCCTGC
 AGACAGCAGTGACGACAGCAGTGACGATTCTGATTCAAGTCTGAGGAAGAGAAAAACCTCCAGCTAAG
 ACGGTCGTCTCCAAGACCCCGCAAAGCAGCTCCAGTGAAGAAGAAAGCAGAAAGCTCTCAGACAGCT
 CGGATTCGACAGTTCTGAGGATGAAGCTCTGCCAAGCCAGTCAGTACAACCAAGAGTCCCAAGCCAGC
 TGTCACTCCGAAGCCATCTGCAGCAAAGGCAGTGACAACCTCTAAGCAACCTGCAGGCAGTAACCCAGAAA
 CCTCAGAGCAGGAAGGCTGACAGCAGCTCCAGCGAGGAGGAAAGCAGCTCCAGCGAGGAGGAGGAGGCT
 CCAAGAAAAGTGCCACAACCCCAAGGCCAAGGTGACTGCTAAAGCAGCACCCGCCAAACAGGCCCTCA
 GGCTGCTGGGACAGCAGCTCTGACTCAGATAGTTCAGCAGTGAAGAGGAGGAGAAGACTCCTAAGCCC
 CCAGCTAAGAAGAAGGCAGCAGGTGGAGCCGTTTCTACACCAGCCCTGGGAAGAAAGCAGAGGCCAAGA
 GCAGCAGCAGCAGCAGCAGCAGCTCCGAAGATTCAGTGAAGAGGAGAAAAAAGAAAGCCAAAGC
 TACTACCCCTAAAATACAGGCAAGCAAGGCCAATGGCACTCCAGTCTCTGAATGGAAAAGCAGCCAAG
 GAAAGTGAGGAGGAGAGGAGGAGGAAGAAACAGAAGAGAAGAAAAGGCAGCTGGGACCAAGCCAGGTT
 CAGGCAAAAAACGAAGCAGAATGAGACCGCAGATGAAGCAACAACCTCCTCAAGCTAAGAAAAGTTAAGCT
 CGAGACCCCAATACGTTTCCAAAAAGGAAGAAGGGAGAAAGAAGGGCGTCTTCCCTTTCCGAAGGGTC
 AGGGAGGAGGAGATTGAGGTGGACTCTCGAGTGGCGGACAATTCTTTGATGCCAAGCGAGGTGCAGCTG
 GAGACTGGGGGAGCGAGCCAATCAGGTTCTGAAGTTCACCAAAGGCAAGTCTTCCGGCATGAAAAAAC
 GAAGAAGAAGCGAGGCAGCTACCGGGAGGCTCCATCTCTGTCCAGGTCAATTCGTCAAATTCGACAGC
 GAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR215487 representing NM_053086
 Red=Cloning site Green=Tags(s)


```
MADTGLRRVPSDL YPL VLRFLRDSQLSEVASKFAKATGATQQDANASSLLDIYSFWLNRSTKAPKVKLQ
SNGPVTKKAKKETSSSDSSEDSEDEDKKAQGLPTQKAAAQVKRASVPQHAGKAAAKASESSSSEESSEE
EEEDKKKKPVQKAAKPQAKAVRPPAKKAESSESDSDSDSDSSSEETPQTQKPKAAVAAKAQTAEAKPG
TPAKAQPKVANGKAAASSSSSSSDDEEKKAAAPPKKTVPKKQVVAKAPVKVAAAPTQKSSSSEDSSS
EEEEGQRQPMKKKAGPYSSVPPPSVPLPKKSPGTQAPKAAAQTQPADSSDDSSDDSDSSSEEEKPPAK
TVVSKTPAKAAPVKKKAESSSDSDSDSSEDEAPAKPVSTTKSPKPAVTPKPSAAKAVTTPKQAPAGSNQK
PQSRKADSSSSEEESSSSEEEEASKKSATTPKAKVTAKAAPAKQAPQAAGDSSSDSDSSSSEEEETPKP
PAKKAAGGAVSTPAPGKAEAKSSSSSSSSSEDSSEEEKKKPKATTPKIQASKANGTPASLNGKAAK
ESEEEEEEEEEETEEKKAAGTKPGSGKRRKQNETADEATTPQAKVKLETPNTPFKRKKGERRASSPFRRV
REEEIEVDSRVADNSFDAKRGAGDWGERANQVLKFTKGKSFRHEKTKKRRGSRGGSISVQVNSVKFDS
E
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



ACCN: NM_053086

ORF Size: 2103 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_053086.3](#), [NP_444316.2](#)

RefSeq Size: 3649 bp

RefSeq ORF: 2106 bp

Locus ID: 70769


UniProt ID: [E9Q5C9](#)

Cytogenetics: 19 C3

MW: 74 kDa

Gene Summary: Nucleolar protein that acts as a regulator of RNA polymerase I by connecting RNA polymerase I with enzymes responsible for ribosomal processing and modification (By similarity). Required for neural crest specification: following monoubiquitination by the BCR(KBTBD8) complex, associates with TCOF1 and acts as a platform to connect RNA polymerase I with enzymes responsible for ribosomal processing and modification, leading to remodel the translational program of differentiating cells in favor of neural crest specification (By similarity). Involved in nucleologenesis, possibly by playing a role in the maintenance of the fundamental structure of the fibrillar center and dense fibrillar component in the nucleolus (PubMed:11424213). It has intrinsic GTPase and ATPase activities (By similarity). [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR215487