

Product datasheet for **MR21524**

Olfr1386 (NM_001011741) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Olfr1386 (NM_001011741) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Olfr1386
Synonyms:	MOR256-50
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR21524 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGATAGTTTCAATGCCACGTTAGAAGAGAGATTCTTTTTGGTGGGCTTTTCAGATTGGCCTCAGCTGG
AACTCATCCTTTTTGTTTTTATTCAATTGTCTACTCCCTAACTATCTTTGGCAACTACCATCATTGC
TCTCTCCGAATAGACCTTCGACTGCACACCCCATGACTTCTCCTCTCAAATCTCTCCTTTCTGGAC
CTCTGCTTACCACCAGCACCGTGCCCGAGCTCCTCATCAACCTCTATGGACAAGACAGGACCATCAGT
ATGGAGGGTGTGTGGCCAGCTGTCATATATCTTGCTCTGGGATCCACTGAGTGTGTGCTCCTGGTAGT
GATGGCTTTTGACCGCTATGCTGCTGTGTGCTGCCACTGCACTACACAACCATCATGCACCCCTTCTC
TGTCCAGGCATTGGCTCTTGCTTCATGGGTAGGAGGATTCTTGAACCTCTGATTCAAACAGGCCTCATGA
TGACCATGCCACTCTGTGGACATCGACTGAATCACTTCTCTGTGAGATGCCTGTCTTCCCAAGTTGGC
CTGCCAGGACACAGGAGGAACAGAGGCCAAGATGTTTGTGGCCAGAGCCATAATCTTGGTCTTCCCTGCA
ACATTAATTCTTGGCTCCTATGGACACATTGCTAAAGCAGTACTAAAGGTCAAGTCAACAGCTGGACGCA
GAAAAGCTTTTGGGACATGTGGTCCCACCTCCTGGTGGTTTTCCCTGTTTTATGGCTCAGCCATTTACAC
ATACTTGCAACCAAGAGCAGTTATTCTGAAAGCAAGGGGAAGTTTGTGGCCCTTTTTATACTATTGTC
ACCCCATGCTCAACCTCTGATCTATACCCTGAGGAACAAGGATGTGAAGGGGGCTCTGTGAAGGTGC
TAGGGAGAGGCACAGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001011741.2](#), [NP_001011741.1](#)

RefSeq Size: 1128 bp

RefSeq ORF: 930 bp


Locus ID: 257888

Cytogenetics: 11 B1.2

MW: 34.4 kDa

Gene Summary: Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for MR215524