

Product datasheet for **MR215568**

Olfr915 (NM_146785) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Olfr915 (NM_146785) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Olfr915
Synonyms:	MOR167-4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR215568 representing NM_146785 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**C

ATGGGTTCGGAATGGCTCTTCAGTGACTGAGTTCATTCTGGTGGGATTAACCAAGGACCCTGATCTCC
 AATGCCCTCTTTTCGTCCTCTTTCTAGTGATGTATGTTCTCACTGTATTGGGAAATCTGAGTTTGATTAC
 CTTAATTGACTGAATTCTCACCTTCATACCCCAATGACTTTTCTCTTTAATTGCTACTAGTTGAC
 TTCTGGTACTCTCTGTGTTACACCAAAAATGCTAATGAGCTTTATATCAGAGAAGAACATTATTTCT
 ATAGAGGATGCATGACCAACTTTCTTTTCTGCTTTTCTGAATATTATGTGTTGACTGT
 AATGGCCTATGATCGCTATGTGGCCATATGTAATCCACTCTGTACAATATTGTCATGTCTCCTAAATTA
 TGTTTGAACCTCATACTTGGCTCATATATAATGGCATTCTAGTGCCATGGCTCACACGGGATGCATGC
 TGAGACTGACCTTCTGTGATGCAACACCATCAATCACTACTTCTGTGACATCCTCCCTTTGCTTCAGCT
 CTCCTGCACCAGCATCTATGTCAACGAGCTGGAGGTTTTGTGTAATAGGCATCAACATTATTGTGCC
 ACCATCACCATCTTTATCTTTATGTTTTATCATCTCAAGTATTTTCGTATCAGCTCCAAGGAGGGCA
 GGTCCAAAGCCTTCAGCACCTGCAGTCCACATAATTGCAGTTTCTCTGTTCTTTGGATCAGGTGCATT
 TATGTATCTCAAACATCTTCTGCTGAGTCAATGAATGAAAGTAAAGTCTCGTCTGTCTTTATACCAAT
 ACAGTTCTCTACTGAATCCCTTAATCTACAGCTTGAGGAACAAAGATGTTAAATGGCTTGATAAAGA
 CCCTGAGTAAGAGAAAGTGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >MR215568 representing NM_146785
 Red=Cloning site Green=Tags(s)


MGSENGSSVTEFILVGLTKDPLQCPLFVLFLVMYVLTVLGNLSLITLIVLNSHLHTPMYFFLFNLSLVD
 FWYSSVFTPKMLMSFISEKNIISYRGCMQLFFFSFFCISEYYVLTVMAYDRYVAICNPLLYNIVMSPKL
 CLNLILGSYIMAFSSAMAHTGCMRLRTFCDANTINHYFCDILPLLQLSCTSIYVNELEVFFVIGINIIVP
 TITIFISYGFIISSIFRISSKEGRSKAFSTCSSHIIAVSLFFGSGAFMYLKPSSAESMNESKVSSVFYTN
 TVPLLNPLIYSLRNKDKIALIKTLSKRKC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9084_g06.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_146785

ORF Size: 930 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_146785.2](#), [NP_666996.2](#)

RefSeq Size: 933 bp

RefSeq ORF: 933 bp


Locus ID: 258781

Cytogenetics: 9 A5.1

MW: 35 kDa

Gene Summary: Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. This olfactory receptor gene appears to represent a strain-specific polymorphic pseudogene in mouse, where some strains, including the C57BL/6 reference genome strain, have a non-functional allele, while other strains have a protein-coding allele. [provided by RefSeq, Jan 2019]

Product images:



Circular map for MR215568