

Product datasheet for **MR215572A1V**

Mouse Olfr159 (NM_019476) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Olfr159 (NM_019476) AAV Particle
Tag:	Myc-DDK
Symbol:	Olfr159
Synonyms:	mOR37e; MOR262-5; Olfr37e
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR215572 representing NM_019476 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGAGAGATCCAACAAGACCACCCCTGTGTCCAGCTTCATCCTCCTGGGTCTCTCCGCCACCCAAAGC
TGGAGAAAACCTTCTTCGTGCTCATCCTGCTGATGTACCTGGTGATCCTGCTGGGCAATGTCGTCCTCAT
CCTGGTGAGCATCCTCGACTCCACCTGCACACGCCCATGTACTTCTCCTGGGGAACCTCTCCTTCCTG
GACATCTGCTACACCACCTCCTCCGTCCCCCTCATTCTGGACAGCTTCCTGACCCCAAGAACCATCTC
CCTTCTCGGGCTGTGCCGTGCAGATGTTTCTCCTTCGCCATGGGAGCCACAGAGTGTGTCTCCTGGG
CATGATGGCGTTTGATCGTTATGTGGCCATCTGCAACCCCTTAGATATCCTGTGGTCATGAGCAAGGCT
GCCTATGTGCCCATGGCCGCGGCTCCTGGGTATCTGGAAGCATCACTGCCACAGTGCAGATATCTTTAG
CAATGACACTGCCCTTCTGTGGGACAATGTCATCAATCACTTTACATGTGAGATCCTGGCAGTCCTGAA
ACTGGCCTGTGCTGACATCTCCATCAATGTCATCAGCATGGCAGTGGCTAATGCAATGTTTTGGGGGTC
CCTGTCCTCTTCATCTTTGTCTCCTACATCTTCATACTCTCCACCATCCTGAGGATCCCCTCTGCCGAGG
GGAGGAAGAAGGCCTTCCACCTGCTCTGCCACCTCACCCTGGTACTTGTCTTCTATGGGACCATCCT
GTTTCATGTACGGAAGCCCAATCAAGGACCATTTGGGGGAGACAAACAGGACCTTGACAGACAAGCTC
ATCTCCCTCTTCTATGGAGTGGTGACCCCATGCTGAACCCCATCATCTACAGCCTGAGGAACAAGGATG
TGAAGGCCGCTGTGACAAACCTGGTGGCCAGAAACACTTCAAGTGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR215572 representing NM_019476 Red=Cloning site Green=Tags(s) MERSNKTTVPVSSFILLGLSAHPKLEKTFVLLILLMYLVILLGNVVLILVSILDSHLHTPMYFFLGNSFL DICYTTSSVPLILDSFLTPTKTIISFGCAVQMFLSFAMGATECVLLGMMAFDRYVAICNPLRYPVVMKA AYVPMAGSWVSGSITATVQISLAMTLPCGDNVINHFTCEILAVLKLACADISINVISMAVANAMFLGV PVLFI FVS YIFILSTILRIPSAEGRKKAFSTCSAHLTVVLV FYGTILFMYGKPKSKDPLGADKQDLADKL ISLFYGVVTPMLNP I IYSLRNKDVKAAVTNLVGQKHFKW TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_019476
ORF Size:	957 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_019476.1, NP_062349.1</u>
RefSeq Size:	960 bp
RefSeq ORF:	960 bp
Locus ID:	29849
Cytogenetics:	4 A5
MW:	35 kDa