

Product datasheet for **MR215590L3V**

Map3k10 (NM_001081292) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Map3k10 (NM_001081292) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Map3k10
Synonyms:	BC028668; BC046514; Mlk2; MST
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001081292
ORF Size:	2826 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215590).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001081292.1 , NP_001074761.1
RefSeq Size:	3690 bp
RefSeq ORF:	2829 bp
Locus ID:	269881
Cytogenetics:	7 A3
Gene Summary:	Activates the JUN N-terminal pathway.[UniProtKB/Swiss-Prot Function]



[View online »](#)