

Product datasheet for **MR215601A1V**

Mouse Nrp (NM_001013372) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Nrp (NM_001013372) AAV Particle
Tag: Myc-DDK
Symbol: Nrp
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR215601 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAATCGAAACCCTGGAGTCGTGACCCCGGAAGAACCTGCCAGAGCCGGAATTCGAGTTCTGCTTCCG
 GGCCAAACTGTTGGCAGCCTCGAGATGGGAAGATGGCGGCTGCTGTGGCTTCATTAGCCACGCTGGCTG
 CAGAGCCAGAGAGGATGCTTTCCGGAAGCTTTTCCGCTTCTACCGGCAGAGCCGGCCGGGACAGCGGA
 CCTGGGAGCCGTCATCGACTTCTCAGAGGCGCACTTGGCTCGGAGCCCGAAGCCCGCGTGTCCCAAGTA
 GGAAAGGAGGAGTAGTGTGTGCCAGCCTAGCGGCCGACTGGGCCACCCGAGACTGGGCCGCTCCGGGCC
 GGCTTTGGAGGAAGCCCTGCTGGGCCTGTCCAGTGAGCTGTAATGTCGAGCGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR215601 protein sequence
 Red=Cloning site Green=Tags(s)

MNRNPGVVTPEEPARAGISSASGPNCWQPRDGEDGGCCGFISHAGCRAQRGCFPEAFPLPAEPAGDSG
 PGSRHRLLRGALGSEPEARAPGRKGGVVCASLAADWATRDWAASGPALEGSPCWACPVSCNVER

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_001013372
ORF Size: 405 bp


[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001013372.2, NP_001013390.1</u>
RefSeq Size:	1170 bp
RefSeq ORF:	408 bp
Locus ID:	654309
Cytogenetics:	12 D2
MW:	13.9 kDa