

## Product datasheet for **MR215635L3V**

### **Erlin1 (NM\_001164359) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Erlin1 (NM_001164359) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Erlin1
Synonyms:	2810439N09Rik; C80197; Keo4; Spfh1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001164359
ORF Size:	1044 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215635).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001164359.1</a> , <a href="#">NP_001157831.1</a>
RefSeq Size:	3245 bp
RefSeq ORF:	1047 bp
Locus ID:	226144
UniProt ID:	<a href="#">Q91X78</a>
Cytogenetics:	19 C3



[View online »](#)

**Gene Summary:**

Component of the ERLIN1/ERLIN2 complex which mediates the endoplasmic reticulum-associated degradation (ERAD) of inositol 1,4,5-trisphosphate receptors (IP3Rs). Involved in regulation of cellular cholesterol homeostasis by regulation the SREBP signaling pathway. Binds cholesterol and may promote ER retention of the SCAP-SREBF complex (By similarity). [UniProtKB/Swiss-Prot Function]