

## Product datasheet for **MR215653L3V**

### **Ppwd1 (NM\_172807) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Ppwd1 (NM_172807) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ppwd1
Synonyms:	4632422M10Rik; A330090G21Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_172807
ORF Size:	1938 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215653).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_172807.4</a> , <a href="#">NP_766395.2</a>
RefSeq Size:	2367 bp
RefSeq ORF:	1941 bp
Locus ID:	238831
UniProt ID:	<a href="#">Q8CEC6</a>
Cytogenetics:	13 D1
Gene Summary:	PPIase that catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides and may therefore assist protein folding. May be involved in pre-mRNA splicing.[UniProtKB/Swiss-Prot Function]



[View online »](#)