

## Product datasheet for **MR215662L3V**

### **Snx10 (NM\_028035) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Snx10 (NM_028035) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Snx10
Synonyms:	2410004M09Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_028035
ORF Size:	606 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215662).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_028035.4</a> , <a href="#">NP_082311.3</a>
RefSeq Size:	2567 bp
RefSeq ORF:	606 bp
Locus ID:	71982
UniProt ID:	<a href="#">Q9CWT3</a>
Cytogenetics:	6 B3



[View online »](#)

**Gene Summary:**

Probable phosphoinositide-binding protein involved in protein sorting and membrane trafficking in endosomes. Plays a role in cilium biogenesis through regulation of the transport and the localization of proteins to the cilium. Required for the localization to the cilium of V-ATPase subunit ATP6V1D and ATP6V0D1, and RAB8A. Involved in osteoclast differentiation and therefore bone resorption.[UniProtKB/Swiss-Prot Function]