

## Product datasheet for **MR215677**

### Lcmt2 (NM\_177846) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Lcmt2 (NM_177846) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Lcmt2
Synonyms:	D330024M17; Tyw4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide  
Sequence:

>MR215677 representing NM\_177846  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGGCCCGAGGCCGCCAGCGCCGGGGAACGGTCCAGAGCACCAACGACAGCAGCTCCCTCAGCA  
 AACGATCGCTTGCCGCACACGGATACGTGCGCGACCCCTTTGCGGCGCTGCTGGTCCCGGACCCGTGCG  
 GCGCACGCCGCTTATCCACCGCGCTATTACGTGCGCGCGCGCCGTGCGTCACTGCGTGCGCCCTTC  
 CTCGAACACTAGCGCGCTCCCGTCCCGGACCAGAGCCAGATCCTGTCTTGGGCTCAGGCTCCGACT  
 CGCTTTATTTTCGTCTGAAAGCGCGGGCTCCTGGCCGAGCTGCCGTCTGGGAGTGGACTTCCCGGA  
 CGTGTCTCGGCTCAAGGCGGAGAGGATCGAGGAGACCCCGGAGCTGCGTGCGCAGACCGTCTTTCAAG  
 ATCGGGGACTCGGCGTACTCTGTGCTTCGAGAGCGCGGACTACCGCATCCTGGGCGGGACCTGCGGG  
 AGCTCCAGCGATTAGGAGAGGCCCTGGATGGCGCCGGCTGGACGCCACCTCTCCACGCTGCTCCTGGC  
 CGAGGCGGTGCTGACCTACTTGGAGCCCTCCAGTGCCACAGCCCTCATTGCTGGGCCGCCAGCGTTTC  
 CCCGACGCCCTTTTCGTGATCTATGAGCAGATGCAGCCGGGGATGCCTTTGGGCAGATCATGCTGCAGC  
 ACTTTCAGCGGTTGCATTACCCTTGCATGGCCTGGAGCTGTTTCCCGTGTGAAGGCCAGCGCCAGCG  
 CTTCTTCAAGCTGGCTGGACTGCCTGCAGCGCCTGGACCTGAACGAGTTCTATCGCCGCTTCTCTCG  
 GCAGAAGAGCGCCAGCGGTTGAGACGCTTGAACCCTTGTGAATATGAGGAGTGGCATCTGAAGTGT  
 CACATTATTTTATCCTGGCAGCATCTAGGGGAGACATTCTGTGAGAACTCCAGTGTGTTGAGCCCTCAGA  
 GGCTTCTTTTTCAGATAGATCCTGCCTCGCCTTCAGGGTTTCTTTCTGCCAGAGTAGTCACTAGTGACCAC  
 CAGCACTCAAGCCTGAAGCGATATGGCCACGCCCTGCGCCTCTGAGCCCTGGCGTTATTTTCAGTGCAG  
 GAGGCTTTGGAGAACAAGAGGGACGACTGCCGAGTGAGCAGGTTTTCAGTGCTCTCCAGATCCTGTGA  
 CTCGGAATGGGAAGGCTGCCAAATAAGTACCTTGGGGACTGAAGGCCAGTGGGATGGACGCTTTATCAC  
 ACCATGACAAGGCTTTCTGATACTCGGGTCTGTTTCTGGAGGAAGACTCTCCCAAGTAAGTCCAGCCT  
 CTGGGGCTCTCCAACCTGATTTATACAAAAGTAAGGATAATTGCTCTGAGGGCCAAAATGTAACAGTAAC  
 GAAGGCTGCCTTGAAGAAGGTTTCAAGTGTGTCATGTTGGCGGCATTCAACAACGGAAGTATACTATCAG  
 AATCAAAGATACTTATTTGTGTATGGTGGCCGAAGTGTGGCTGAACCTGTACTAAGTACTGTCGTTTTTC  
 TACATGTAGAGACAATGGCTTGGGTCAGAATCCCAGTACAAGGCGCATCACCCGAAGGTCGGCATTCCCA  
 TAGCGCCTGCAGTTGGCAAGGGGAGCACTTATCGCTGGAGGTCTTGGGGCTTCTGAGGAACTGTTGAGC  
 TCTGTACTCTTTCTCAAACAGTGTCTCTGGATTCTCTGGGAATCAATAGATATCCAGCCCTCCATTA  
 CCCAAGTACTCTCACACGGCTCATGTGTTAATGGAAAGCTGCTTCTGGTTGGAGGGTCTGGATTCA  
 TTCTCTCAGTTCCTGGAGTACTGTTATCTGTTTACTACAGGGTTGAGCTCTGAATATCAGATTGAC  
 ACAGCGTCTGTGCCATGGCCATTAATGCTGCACAACCACAGCAGCGCCCTCCTTCTGAAGAGCAGCAGC  
 TCTGTATTGGAGGTGGAGGAACTGCTTTTCTTTGGTACTTATTTCAACCCACACAGTGGCATT  
 AGACCTTCTTCTTACGTTTACGGCAA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR215677 representing NM\_177846  
Red=Cloning site Green=Tags(s)


MGPRGRQRRAGTVQSTNDSSSLSKRSLAAHGYVRDPFAALLVPGPVRRTPLIHRGYYVRARAVRHCVRAF  
 LELTSALPSRTRAQILSLGSGSDSLYFRLKAAGLLARAAVWEVDFPDVSRKKAERIEETPELRAQTGPFK  
 IGDSSSLCFESADYRILGADLRELQRLGEALDGAGLDATSPTLLLAEAVLTYLEPSSATALIAWAAQRF  
 PDALFVIYEQMQPGDAFGQIMLQHFQRLHSPHLGLELFPVVKARQRFLQAGWTACSLDLNEFYRRLS  
 AEERQRVETLEPFDEYEEWHLKCSHYFILAASRGDILSETPVFEPSEASFQIDPASPGLSARVVTSDH  
 QHSSLKRYGHASALLSPGVIFSAGGFGEQEGRHCRVSRFHVLSRSCDSEWEGCQISTLGTEGQWDGRLYH  
 TMTRLSOTRVLVGGRLSPVSPASGALQLDLKSKDNCSEGNVTVTKAALEEGSVLSCWRHSTTEVYYQ  
 NQRYLFVYGGRSVAEPLSDCRFLHVETMAWVRIPVQGASPEGRHSHSACSQQGGALIAGGLGASELLS  
 SVLFLKPVSSGFLWESIDIQPSITPRYSHTAHVFNKLLLVGGVVIHSSSVPGVTVICLTTGLSSEYQID  
 TASVPWPLMLHNHSSALLPEEQQLLLIGGGGNCFSFGTYFNPHTVALDLSSLRSGQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mm9100\\_a07.zip](https://cdn.origene.com/chromatograms/mm9100_a07.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_177846


**ORF Size:** 2058 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u>NM_177846.3, NP_808514.2</u>
<b>RefSeq Size:</b>	3407 bp
<b>RefSeq ORF:</b>	2061 bp
<b>Locus ID:</b>	329504
<b>UniProt ID:</b>	<u>Q8BYR1</u>
<b>Cytogenetics:</b>	2 E5
<b>MW:</b>	75.3 kDa
<b>Gene Summary:</b>	Probable S-adenosyl-L-methionine-dependent methyltransferase that acts as a component of the wybutosine biosynthesis pathway. Wybutosine is a hyper modified guanosine with a tricyclic base found at the 3'-position adjacent to the anticodon of eukaryotic phenylalanine tRNA (By similarity). May methylate the carboxyl group of leucine residues to form alpha-leucine ester residues.[UniProtKB/Swiss-Prot Function]

### Product images:



Circular map for MR215677