

Product datasheet for **MR215677L4V**

Lcmt2 (NM_177846) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lcmt2 (NM_177846) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Lcmt2
Synonyms:	D330024M17; Tyw4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_177846
ORF Size:	2058 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215677).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_177846.3 , NP_808514.2
RefSeq Size:	3407 bp
RefSeq ORF:	2061 bp
Locus ID:	329504
UniProt ID:	Q8BYR1
Cytogenetics:	2 E5



[View online »](#)

Gene Summary:

Probable S-adenosyl-L-methionine-dependent methyltransferase that acts as a component of the wybutosine biosynthesis pathway. Wybutosine is a hyper modified guanosine with a tricyclic base found at the 3'-position adjacent to the anticodon of eukaryotic phenylalanine tRNA (By similarity). May methylate the carboxyl group of leucine residues to form alpha-leucine ester residues.[UniProtKB/Swiss-Prot Function]