

## Product datasheet for MR215702

### Ndufa3 (NM\_025348) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Ndufa3 (NM\_025348) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Ndufa3  
**Synonyms:** 1010001M12Rik; 1700022J01Rik  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >MR215702 representing NM\_025348  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCCGGGAGAATCTCTGCCTTCTCAAGAATGCCTGGGCGAAGGAGCCGGTGCTGGTGGTGCCTTCT  
CTGTCTGGGGCTCGCTATAATTATGCCCATGATTAGCCCTACACCAAGTATGCTAGCATGATCAACAA  
GGCCACACCTACAACCTACCCAGTGCCTGTGAGAGATGACGGGAACATGCCTGATGTGCCAGCCACCT  
CAGGATCCTTTGGGTCCAAGCTTGGACTGGCTGAAGAACCTG

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR215702 representing NM\_025348  
Red=Cloning site Green=Tags(s)

MAGRISAFLNKAWAKEPVLVVSFSVWGLAIIMPDISPYTKYASMINKATPYNYPVPRDDGNMPDVPSPH  
QDPLGPSLDWLKNL

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

**Chromatograms:** [https://cdn.origene.com/chromatograms/mm9026\\_d02.zip](https://cdn.origene.com/chromatograms/mm9026_d02.zip)

**Restriction Sites:** Sgfl-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_025348

**ORF Size:** 252 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_025348.2](#), [NP\\_079624.1](#)

**RefSeq Size:** 349 bp

**RefSeq ORF:** 255 bp

**Locus ID:** 66091


UniProt ID: [Q9CQ91](#)

Cytogenetics: 7 A1

MW: 9.8 kDa

Gene Summary: Accessory subunit of the mitochondrial membrane respiratory chain NADH dehydrogenase (Complex I), that is believed not to be involved in catalysis. Complex I functions in the transfer of electrons from NADH to the respiratory chain. The immediate electron acceptor for the enzyme is believed to be ubiquinone.[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for MR215702