

Product datasheet for **MR215722L3V**

Nap1l2 (NM_008671) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	Nap1l2
Synonyms:	Bpx
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_008671
ORF Size:	1380 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR215722).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_008671.2
RefSeq Size:	2606 bp
RefSeq ORF:	1383 bp
Locus ID:	17954
UniProt ID:	P51860
Cytogenetics:	X 45.97 cM



Gene Summary: Acidic protein which may be involved in interactions with other proteins or DNA.
[UniProtKB/Swiss-Prot Function]