

Product datasheet for **MR215724L3V**

Calr3 (NM_028500) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Calr3 (NM_028500) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Calr3
Synonyms:	1700031L01Rik; 6330586I20Rik; Crt2; cspn
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_028500
ORF Size:	1140 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215724).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_028500.3 , NP_082776.2
RefSeq Size:	1670 bp
RefSeq ORF:	1143 bp
Locus ID:	73316
UniProt ID:	Q9D9Q6
Cytogenetics:	8 B3.3
Gene Summary:	CALR3 capacity for calcium-binding may be absent or much lower than that of CALR (By similarity). During spermatogenesis, may act as a lectin-independent chaperone for specific client proteins such as ADAM3. Required for sperm fertility.[UniProtKB/Swiss-Prot Function]



[View online »](#)