

Product datasheet for **MR215749L3V**

Chsy1 (NM_001081163) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Chsy1 (NM_001081163) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Chsy1
Synonyms:	mKIAA0990
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001081163
ORF Size:	2400 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215749).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001081163.1
RefSeq Size:	4176 bp
RefSeq ORF:	2403 bp
Locus ID:	269941
UniProt ID:	Q6ZQ11
Cytogenetics:	7 C



[View online »](#)

Gene Summary:

Has both beta-1,3-glucuronic acid and beta-1,4-N-acetylgalactosamine transferase activity. Transfers glucuronic acid (GlcUA) from UDP-GlcUA and N-acetylgalactosamine (GalNAc) from UDP-GalNAc to the non-reducing end of the elongating chondroitin polymer.
[UniProtKB/Swiss-Prot Function]