

Product datasheet for **MR215755L3V**

Egfl7 (NM_198725) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Egfl7 (NM_198725) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Egfl7
Synonyms:	VE-statin; Zneu1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_198725
ORF Size:	795 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215755).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_198725.2 , NP_942018.2
RefSeq Size:	1365 bp
RefSeq ORF:	798 bp
Locus ID:	353156
UniProt ID:	Q9QXT5
Cytogenetics:	2 A3
Gene Summary:	Regulates vascular tubulogenesis in vivo. Inhibits platelet-derived growth factor (PDGF)-BB-induced smooth muscle cell migration and promotes endothelial cell adhesion to the extracellular matrix and angiogenesis.[UniProtKB/Swiss-Prot Function]



[View online »](#)