

## Product datasheet for **MR215769**

### Lingo2 (NM\_001166001) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Lingo2 (NM_001166001) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Lingo2
Synonyms:	B230217C06Rik; Lern3; Lrrn6c
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>MR215769 ORF sequence  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGCTTCACACGGCTATACCATGCTGGCAGCCATTCTGGGTCTGGCTGTGGTGTACTCTTAATGGGAT  
 CCACCATTGGCTGCTCTGCTCGCTGTGAGTGCTCCGCCAGAACAATCTGTTAGCTGCCACAGAAGACG  
 ATTGCTCGCGATCCCAGAAGGCATTCCATTGAGACCAAATCTTGGACCTGAGCAAAAATCGACTAAAG  
 AGCATAAACCTGAAGAGTTCATCTCATCTCTGTTGGAGGAGATAGACTTGAGCGACAACATTATTG  
 CCAATGTGGAGCCTGGGGCATTTAAACATCTTTAACCTGCGTTCCTCCGCCTAAAAGGCAATCGCCT  
 TAAGTTGGTCCCTTAGGAGTATCACAGGACTGTCCAACCTCACCAAGCTTGACATTAGTGAGAATAAG  
 ATTGTCATTTTGTGGACTACATGTTCCAGGATCTGCATAACCTGAAGTCTCTAGAAGTGGGGGACAATG  
 ATTTAGTGTATATCTCACACAGGGCCTCAGCGGACTACTTAGCTGGAGCAGCTCACCTGGAGAAGTG  
 CAACTTGACAGCAGTACCAACAGAAGCCCTTTCCCATCTCCGAGCCTCATCGCCCTGCATCTGAAGCAT  
 CTCAATATCAACAATATGCCTGTGTATGCCTTTAAAAGATTGTTCCACCTGAAAAACCTAGAGATCGACT  
 ATTTGGCCTTTGTTGGATTGATGCCAGCCAACAGCCTCTATGGTCTCAACCTCACGTCCCTTTCAATCAC  
 CAACACCAACCTGTCCACTGTCCCTTCTCGCCTTTAAACACCTTGATACCTGACCCACCTTAACCTC  
 TCCTACAATCCCATCAGCACTATTGAAGCTGGCATGTTCTCTGACCTGATCCGCCTACAGGAGCTTCATA  
 TAGTGGGGGCCAGCTCCGCACTATTGAGCCTCACTCCTTCCAAGGGCTCCGCTTCTCCGTGTGCTCAA  
 TGTATCTCAGAACCTGCTGGAACATTGGAAGAGAAGCTTTCTCCTCCCTAGGGCTTTGGAGGTCCTG  
 AGCATTAAACAACAACCCACTAGCCTGTGACTGCCGACTCCTCTGGCTCCTGCAGCGACAACCCAACCTGC  
 AGTTTGGGGGCCAGCAGCCCATGTGTGCTGGGCCAGACACCATCCGTGAGAGATCATTTAAGGATTTCCA  
 TAGCACTGCTCTTTCTTTTATTTTACCTGCAAAAAACCAAAAATCCGTGAAAAGAAGTTACAGCATCTC  
 CTCGTGGATGAAGGGCAAACGGTCCAGCTGGAGTGCAACGCTGATGGAGACCCGAGCCCGTGATTTCT  
 GGGTGACACCTCGAAGGCGTTTTATCACACCAAGTCCAACGGAAGGGCCACTGTGTTGGGTGATGGCAC  
 CTTGGAAATCCGTTTTGCCAGGATCAAGACAGTGGGATGATGTTTGCATCGCGAGCAACGCTGCTGGG  
 AACGATACCTTACAGCATCTCTACTGTGAAGGGATTACGTCAGACCGCTTCTTTACGCAAACAGGA  
 CCCCTATGTACATGACTGACTCCAACGACACCGTTTCCAACGGCACTAATGCCAATACTTTCTCCCTGGA  
 CCTTAAAAAATACTGGTATCTACAGCCATGGGCTGTTTACATTCCTGGGAGTGGTTTTATTTGTTTT  
 CTCCTTCTTTTGTGTGGAGCCGAGGGAAGGCAAGCACAAAAACAGATTGACCTTGAGTATGTGCCCC  
 GAAAAACAATGGTGCTGTTGTGGAAGGGGAGGTGGCTGGCCCCAGGAGTTCAACATGAAAATGATC

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:**

>MR215769 protein sequence  
 Red=Cloning site Green=Tags(s)

MLHTAIPCWQPFLLGLAVVLLMGSTIGPCARCECSAQNKSVSCHRRRLLAIPGPIETKILDLSKNRLK  
 SINPEEFISYPLLEEIDLSNIIANVEPGAFNNLFNLRSLRLKGNRLKLVPLGVFTGLSNLTKLDISENK  
 IVILLDYMFDLHNLKSLEVDNDLVYISHRAFSGLLSLEQLTLEKCNLTAVPTEALSHLRSLIALHLKH  
 LNINMPVYAFKRLFHLKNLEIDYWPLDLMPANSLYGLNLTSLITNTNLSTVPFLAFKHLVYLTHLNL  
 SYNPISTIEAGMFSDLIRLQELHIVGAQLRTIEPHSFQGLRFLRVLNVSQNLLETLEENVFSSPRALEVL  
 SINNNPLACDCRLLWLLQRPNLQFGGQPPMCAGPDTIRERSFKDFHSTALSFYFTCKKPKIREKKLQHL  
 LVDEGQTVQLECNADGDPQPVISWVTPRRRFITTKSNGRATVLDGDTLEIRFAQDQDSGMYVCIASNAAG  
 NDFTTASLTVKGFSDRFLYANRTPMYMTSDNDSVSNGTNANTFSLDLKTIILVSTAMGCFTFLGVVLF  
 LLLFVWSRGKHKNSIDLEYVPRKNGAVVEGEVAGPRRFNMKMI

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:**

Sgfl-MluI

Cloning Scheme:



ACCN: NM\_001166001

ORF Size: 1821 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.


Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001166001.1](#), [NP\\_001159473.1](#)

RefSeq Size: 3562 bp  
 RefSeq ORF: 1821 bp  
 Locus ID: 242384  
 UniProt ID: [Q3URE9](#)  
 Cytogenetics: 4 A5  
 MW: 68.1 kDa

**Product images:**



Circular map for MR215769