

Product datasheet for **MR215774L3V**

Alg1 (NM_145362) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Alg1 (NM_145362) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Alg1
Synonyms:	HMAT1; HMT1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_145362
ORF Size:	1446 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215774).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145362.2 , NP_663337.2
RefSeq Size:	1772 bp
RefSeq ORF:	1449 bp
Locus ID:	208211
UniProt ID:	Q921Q3
Cytogenetics:	16 A1
Gene Summary:	Participates in the formation of the lipid-linked precursor oligosaccharide for N-glycosylation. Involved in assembling the dolichol-pyrophosphate-GlcNAc(2)-Man(5) intermediate on the cytoplasmic surface of the ER (By similarity).[UniProtKB/Swiss-Prot Function]


[View online »](#)