

Product datasheet for **MR215787L4V**

Gabrg3 (NM_008074) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gabrg3 (NM_008074) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gabrg3
Synonyms:	B230362M20Rik; Gabrg-3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_008074
ORF Size:	1401 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215787).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_008074.2
RefSeq Size:	1673 bp
RefSeq ORF:	1404 bp
Locus ID:	14407
UniProt ID:	P27681
Cytogenetics:	7 33.48 cM
Gene Summary:	GABA, the major inhibitory neurotransmitter in the vertebrate brain, mediates neuronal inhibition by binding to the GABA/benzodiazepine receptor and opening an integral chloride channel.[UniProtKB/Swiss-Prot Function]



[View online »](#)