

## Product datasheet for **MR215797L3V**

### Stx16 (NM\_001102424) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Stx16 (NM_001102424) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Stx16
Synonyms:	4930401D03; 5430410K23Rik; 6330500A18Rik; A1648908; AW553605; Syn16
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001102424
ORF Size:	927 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215797).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001102424.1</a> , <a href="#">NP_001095894.1</a>
RefSeq Size:	3717 bp
RefSeq ORF:	930 bp
Locus ID:	228960
Cytogenetics:	2 H4
Gene Summary:	SNARE involved in vesicular transport from the late endosomes to the trans-Golgi network. [UniProtKB/Swiss-Prot Function]



[View online »](#)