

Product datasheet for **MR215800A1V**

Mouse Dadl (NM_001113358) AAV Particle

Product data:

Product Type: AAV Particles

Tag: Myc-DDK

Symbol: Dadl

Synonyms: AI323713

Mammalian Cell None

Selection:

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >MR215800 representing NM_001113358
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGTCGGCGTCTGTGGTGTCCGTCACTCCCGTTCTGGAGGAGTACTTGAGCTCCACTCCGCAGCGGC
TGAAGTTGCTGGACGCCTATCTCTTTATATACTGCTGACCGGGCGCTGCAGTTCGGCTACTGTCTCT
CGTGGGCACCTTCCCTTCAACTCTTCTCTCTGCTTCACTCTGTGTGGGCAGCTTCACTCCTAGCG
GTTTGCTGAGAATACAGATCAACCCCAAGCAAGGCGGACTTCAAGGCATCTCTCCGAGCGAGCCT
TTGCTGACTTCTCTTGGCAGCAGATCTGCACCTTGTCGTCATGAACCTCGTTGGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR215800 representing NM_001113358
Red=Cloning site Green=Tags(s)

MSASVSVISRFLEEYLSSTPQRLKLLDAYLLYILLTGALQFGYCLLVGTFFPNSFLSGFISCVGSFIL
VCLRIQINPQNKADFQGISPERAFADFLFASTILHLVVMNFVG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse

Serotype: AAV-2

ACCN: NM_001113358

ORF Size: 339 bp



Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_001113358.1 , NP_001106829.1
RefSeq Size:	713 bp
RefSeq ORF:	342 bp
Locus ID:	13135
UniProt ID:	P61804
Cytogenetics:	14 27.7 cM
MW:	12.5 kDa