

Product datasheet for **MR215801A1V**

Mouse *Snrpf* (NM_027246) AAV Particle

Product data:

Product Type: AAV Particles
Tag: Myc-DDK
Symbol: *Snrpf*
Synonyms: 2010013O18Rik; sm-F; SMF; snRNP-F
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >MR215801 representing NM_027246
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**C

ATGAGTCTACCCCTCAACCCCAAGCCCTTTCTCAATGGGCTGACAGGGAAGCCAGTGATGGTGAACTTA
AGTGGGCGATGGAGTACAAGGGCTACCTGGTCTCCGTGGATGGCTACATGAACATGCAGCTTGCAATAC
AGAAGAGTACATAGATGGGGCGCTGTCTGGACATCTGGCGAAGTTCTAATAAGGTGTAATAATGTCCTC
TACATCAGAGCGTTGAAGAAGAAGAGGAAGATGGGGAATGAGAGAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR215801 representing NM_027246
Red=Cloning site Green=Tags(s)

MSLPLNPKPFLNGLTGKPMVKLKWGMEYKGYLVSDGYMNMQLANTEEYIDGALSGHLGEVLIRCNVL
YIRGVEEEEEDGEMRE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_027246
ORF Size: 258 bp



Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_027246.1, NP_081522.1</u>
RefSeq Size:	869 bp
RefSeq ORF:	261 bp
Locus ID:	69878
UniProt ID:	<u>P62307</u>
Cytogenetics:	10 C2
MW:	10.2 kDa