

Product datasheet for **MR215808L4V**

Shc2 (NM_001024539) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Shc2 (NM_001024539) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Shc2
Synonyms:	6720466E06; SCK; ShcB; Sli
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001024539
ORF Size:	1719 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215808).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001024539.1 , NP_001019710.1
RefSeq Size:	3974 bp
RefSeq ORF:	1722 bp
Locus ID:	216148
UniProt ID:	Q8BMC3
Cytogenetics:	10 39.72 cM
Gene Summary:	Signaling adapter that couples activated growth factor receptors to signaling pathway in neurons. Involved in the signal transduction pathways of neurotrophin-activated Trk receptors in cortical neurons.[UniProtKB/Swiss-Prot Function]



[View online »](#)