

## Product datasheet for **MR215835**

### **Mtmr7 (NM\_001040699) Mouse Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Mtmr7 (NM_001040699) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Mtmr7
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide  
Sequence:

>MR215835 representing NM\_001040699  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC

ATGGAGCACATCCGCACGCCAAGTTGAAAACGTGCGCTTGGTAGACAGAGTGTCTTGTAAAAAGGCAG  
CGCTGGGAACTTTGTATCTGACGGCCACTCACGTCATATTCGTGGAAAATGCACCTGACACAAGAAAAGA  
AACCTGGATTCTGCACAGTCAGATCTCCACCATTGAAAAGCAGGGCACAACAGCTACGGGCTGCCCTCTG  
CTTATTCGCTGCAAGAACTTCCAGATCGTACAACCTTGTCACTCCACAGGAGAGAGACTGCCACGACGTGT  
ACATCTCTCTGATACGCCTTGCACGGCCAGTGAAGTACGAGGAGCTATACTGCTTTTCTTCAACCCCAA  
GCTGGATAAGGAAGAACGGGAGCAGGGCTGTTGCTGTTGACCTTAGTGAAGAATAACAAGCAATGGGC  
CTCCCTGATAATTATTGGCAGCTCAGCGATGTGAATAGAGACTACAGAGTTTGTGACTCTTACCCACGG  
AACTGTACGTCCCAGATCGGCCACGGCACACATCATCGTGGCAGCTCCAAGTTCGGAGCAGACGGCG  
ATTTCTGCCCTTCTACTACTGTAAAGACAGCCATGCCTCCATCTGCAGGAGCAGCCAGCCCCTTTCC  
GGCTTCAGCGCTCGGTGCCTGGAGGATGAGCAGATGCTCCAAGCCATCAGGAAAGCCAATCCTGGAAGTG  
ACTTCATTTATGTCGTGGACACAAGGCCTAAACTTAATGCGATGGCAAAATCGTGTGCAGGGAAAGGCTA  
TGAGAATGAAGACAATTATTCCAATATCAAGTTTCACTTATCGGCATAGAGAACATCCACGTATGAGG  
AACAGTCTTCAGAAAATGCTGGAAGTGTGTGAACCTAAATCTCCCTCTATGAGTGATTTCTGTGGGGAC  
TGGAGAACTCAGGCTGGCTGAGGCACATCAAAGCCATCATGGATGCAGGGATCTTCATTGCAAAGGCAGT  
TTCAGAGGAAGGGCAAGTGTGCTTGTCTACTGTTCTGACGGCTGGGACAGAACGGCCAGGTGTGCTCG  
GTGGCTAGCCTGCTGCTGGACCCCTACTACCGACTCTGAAAGGATTCATGGTATTAATTGAAAAGGACT  
GGATTTCTTCGGACATAAATTTAATCATAGATATGGCAACTTAGATGGTGACCCAAAAGAAATCTCTCC  
AGTAATCGACCAGTTCATTGAGTGTGTCTGGCAGCTAACGGAACAGTTCCTTGCCTTTGAGTTCAAC  
GAGAGTTTCTGACTCACATTCAGCATCATGTTTATTCTGCCAGTTTGGGAACTTCTGTGTAATAGCC  
AGAAGGAGAGGCGAGAATCAAGATTCAAGAAAGAACCTACTCTTGTGGTCCAACCTGTGGAAGAATCG  
AGCTGACTACCTCAATCCTCTGTTTAGAGCTGACCACAGTCAGACCCAGGGGAGCCTTCTCTGCCTACA  
GCTCCATGCAACTTCACTATAAGTTTTGGAATGGGATGTATAACCGCTTTGAAAAGGGTTGCAGCCTC  
GACAGTCTGTTACAGATTACCTCATGGCTGTGAAGGAGGAGTCCCAGCAGCTGGAGGAAGAACTGGAGTC  
CCTAGAAGAAAGGCTGGAAAAAATTCAGAAAGTCCACTTACATGGTACCAAAGTGAAGAGTAAGCAAAGT  
GAGCCCAGCAAACACTCGGGTTCTCTACCTCAGACCACAGCACAGCCAACACCCTCAGGATTACAGTG  
GGAACAGCAAGTCTTCCATCCAGGAGCCCGTCCCAGGGTATGAAGATTCTGCTCTGATTCTAACCCA  
AGACAATTTAAAAAGCTCAGACCCCGACCTCTCAGTGAATAGTGATCAAGAGTCTGGGGTGGAGGATTTG  
AGCTGTGATCTCCAGTGGTGGTGAAGTGCACCCAGCGAAGACAGTGGCAAGGACCGGGACTCGGACG  
AAGCTGTGTTTCTCACTGCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR215835 representing NM\_001040699  
 Red=Cloning site Green=Tags(s)


MEHIRTPKVENVRLVDRVSKKAALGTL YLTATHVIFVENAPDTRKETWILHSQISTIEKQATTATGCPL  
 LIRCKNFQIVQLVIPQERDCHDVYISLIRLARPVKYEEL YCF SFNPKLDKEEREQGWLLVDLSEEYKRMG  
 LPDNYWQLSDVNRDYRVCD SYPTELYVPRSAHAHIVGSSKFRSRRFPALSY YCKDSHASICRSSQPLS  
 GFSARCLEDEQMLQAIRKANPGSDFIYVVDTRPKLNAMANRAAGKGYENEDNYSNIKFQFIGIENIHVMR  
 NSLQKMLEVCELKSPMSDFLWGLENSGWL RHIKAIMDAGIFIAKAVSEEGASVLVHCSDGWDRTAQVCS  
 VASLLLDPPYYRTLKGFVMLIEKDWISFGHKFNHRYGNLDGDPKEISPVIDQFIECVWQLTEQFPCAFEFN  
 ERFLTHIQHHVYSCQFGNFLCNSQKERRELKIQERTYSLWSNLWKNRADYLNPLFRADHSQTQGSLLHPT  
 APCNFTYKFWNGMYNRFKGLQPRQSVTDYLMAYKEESQQL EEELESLEERLEKIQKVHLHGTVKSKQS  
 EPSKHSGFSTSDHSTANTPQDYSGNSKSFSPSRSPSQGDEDSALIL TQDNLKSSDPDL SVNSDQESGVEDL  
 SCRSPSGGEHAPSEDSGKDRDSDEAVFLTA


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



**Plasmid Map:**


**ACCN:** NM\_001040699

**ORF Size:** 1980 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001040699.1](#), [NP\\_001035789.1](#)  
RefSeq Size: 3446 bp  
RefSeq ORF: 1983 bp  
Locus ID: 54384  
UniProt ID: [Q9Z2C9](#)  
Cytogenetics: 8 A4  
MW: 76.1 kDa  
Gene Summary: Phosphatase that specifically dephosphorylates phosphatidylinositol 3-phosphate (PtdIns(3)P) and inositol 1,3-bisphosphate (Ins(1,3)P2).[UniProtKB/Swiss-Prot Function]