

Product datasheet for **MR215839L3V**

Nt5c3b (NM_001102650) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nt5c3b (NM_001102650) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Nt5c3b
Synonyms:	2610037D24Rik; AI841843; C330027I04Rik; Nt5c3l
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001102650
ORF Size:	900 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215839).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001102650.1
RefSeq Size:	1416 bp
RefSeq ORF:	903 bp
Locus ID:	68106
UniProt ID:	Q3UFY7
Cytogenetics:	11 D



[View online »](#)

Gene Summary:

Specifically hydrolyzes 7-methylguanosine monophosphate (m(7)GMP) to 7-methylguanosine and inorganic phosphate. The specific activity for m(7)GMP may protect cells against undesired salvage of m(7)GMP and its incorporation into nucleic acids. Also has weak activity for CMP. UMP and purine nucleotides are poor substrates (By similarity).[UniProtKB/Swiss-Prot Function]