

Product datasheet for MR215862

Otor (NM_020595) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Otor (NM_020595) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Otor
Synonyms: CDRAP; Fdp; MIA; MIAL
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR215862 representing NM_020595
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCAAGGATATTGATTCTTTGCTTGGGGCCCTTGTGGTTCTATGTGCCGGGCATGGTGTATTTATGG
 AATAACTTTCTTCTAAGAAGTTGTGTCCGATGAGGAGTGTGTCTATACTATTTCTCTGGCAAGAGCACA
 GGAAGATTACAATGCCCCAGACTGTAGGTTTCATCGATGTCAAGAAAGGCAGCAGATCTATGTTTACTCC
 AAGCTGGTAACAGAAAACGGAGCTGGAGAGTTTGGGCTGGCAGTGTATGTTGACCACCAGGATGAGA
 TGGGAATTGTAGTTATTTCCCAAGCACTTGGTGAAGGAGCAGCGTGTATACCAGGAGCCACCAAGGA
 GATCCCAACCACGGATATTGACTCTTCTGTGAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR215862 representing NM_020595
 Red=Cloning site Green=Tags(s)

MARILILLGGLVVLCAHGVMFKLSSKKLCADEECVYIISLARAQEDYNAPDCRFIDVKKGQQIYVYS
 KLVTEGAGEFWAGSVYGDHQDEMGIYGFPSNLVKEQRYVQEATKEIPTDIDFFCE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9006_d02.zip

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:



ACCN: NM_020595

ORF Size: 384 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_020595.3](#)

RefSeq Size: 929 bp


RefSeq ORF: 387 bp

Locus ID: 57329

UniProt ID: [Q9JIE3](#)

Cytogenetics: 2 G1
 MW: 14.8 kDa

Product images:



Circular map for MR215862