

## Product datasheet for **MR215868L4V**

### Tlk2 (NM\_011903) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Tlk2 (NM_011903) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tlk2
Synonyms:	4933403M19Rik; PKUalpha; Tlk
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_011903
ORF Size:	2154 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215868).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_011903.3</a> , <a href="#">NP_036033.2</a>
RefSeq Size:	5379 bp
RefSeq ORF:	2157 bp
Locus ID:	24086
Cytogenetics:	11 E1



[View online »](#)

**Gene Summary:**

Serine/threonine-protein kinase involved in the process of chromatin assembly and probably also DNA replication, transcription, repair, and chromosome segregation. Phosphorylates the chromatin assembly factors ASF1A AND ASF1B. Phosphorylation of ASF1A prevents its proteasome-mediated degradation, thereby enhancing chromatin assembly (By similarity). Negative regulator of amino acid starvation-induced autophagy (By similarity). [UniProtKB/Swiss-Prot Function]