

Product datasheet for **MR215880L3V**

Dcaf6 (NM_028759) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Dcaf6 (NM_028759) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Dcaf6
Synonyms:	1200006M05Rik; lqwd1; PC326
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_028759
ORF Size:	2628 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215880).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_028759.1 , NP_083035.1
RefSeq Size:	3306 bp
RefSeq ORF:	2631 bp
Locus ID:	74106
UniProt ID:	Q9DC22
Cytogenetics:	1 H2.2- H2.3
Gene Summary:	Ligand-dependent coactivator of nuclear receptors. Enhance transcriptional activity of the nuclear receptors NR3C1 and AR. May function as a substrate receptor for CUL4-DDB1 E3 ubiquitin-protein ligase complex (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)