

Product datasheet for **MR215890L3V**

Serbp1 (NM_025814) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Serbp1 (NM_025814) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Serbp1
Synonyms:	1200009K13Rik; 9330147J08Rik; AL022786; Pairbp1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_025814
ORF Size:	1221 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215890).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025814.2 , NP_080090.2
RefSeq Size:	6685 bp
RefSeq ORF:	1224 bp
Locus ID:	66870
UniProt ID:	Q9CY58
Cytogenetics:	6 C1
Gene Summary:	May play a role in the regulation of mRNA stability. Binds to the 3'-most 134 nt of the SERPINE1/PAI1 mRNA, a region which confers cyclic nucleotide regulation of message decay. Seems to play a role in PML-nuclear bodies formation.[UniProtKB/Swiss-Prot Function]



[View online »](#)