

## Product datasheet for **MR215910**

### **Btbd11 (NM\_001017525) Mouse Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Btbd11 (NM_001017525) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Btbd11
Synonyms:	6330404E16Rik
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide  
Sequence:

>MR215910 representing NM\_001017525  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGAGGAAACCGAGGCCAAGAGAGTCTAGAGAGCCTGACCCGGGGTCACCAGTGGGGTCAGGATGTCTCC  
AAAGGCTCTCGCACTCAGGAGTTTCCAGGCCCTCCTCGTGCTGCCCCCTCATGGAGTGGATCCGAGT  
GGCGTGCGCATGCCGGCCACCGCCGAGCTTTTCTATGGATAGCGACGATGTCCGCCAAGCCGCCCGG  
CTGCTGCTCCGGGAGTGGACTGCGAGCCGCGCCAGCTCAAGGCGGACGACTGCTTTTGTCTTCTCGGA  
AGCTGGACGCCGTGGCCATCGAAGCCAAGTTCAAGCAGGACTTGGGCTTCCGGATGCTGAATTGTGGACG  
GACCGACTGGTGAAGCAGGCAGTGTCTTTGCTGGGGCCGGATGGAATCAACACCATGAGCGAACAGGGT  
ATGACCCCTTGATGTACGCCTGCGTCCGTGGGACGAGGCAATGGTGCAGATGCTGCTGGATGCCGGAG  
CTGACCTGAATGTGGAGTTGTCACTCCTCATAAATATCCATCCGTCCACCCCGAGACCCGCCATTG  
GACGGCTCTGACTTTTGTGTGTGCATGGACATATTCCTGTAGTTCAGCTCCTCCTGGATGCTGGGGCC  
AAGGTGGAAGGCTCTGTGGAGCATGGCGAGGAGAACTATTCAGAACTCCCTCCAGCTTGACGCCGCTG  
TTGAAATTTTGGAGTGGTTAGTTTGTGTTGGAGCGTGGTGGGATCCCTGATAGGAACCATGTACAG  
GAATGGAATTTCTACAACCCCCAGGGTATGAACTCTTTCAGCCAGGCCGACGCCATGGACACAGG  
AATGTGTTCCGCAAACTGCTCGCTCAGCCAGAGAAGGAGAAGAGCGATATCCTGTCCCTGGAGGAGATTC  
TGGCTGAGGGGACTGACCTGGCAGAAACAGCCCCGCCCTTGTGTGCCAGCCGCAACAGCAAGGCCAA  
ACTGAGGGCCCTGAGGGAGGCCATGTATCACAGCGCTGAGCACGGCTACGTGGATGTACAAATTGATATC  
AGGAGTATAGGTGTCCCATGGACTCTACACACATGGCTGGAGTCACTTCGGATTGCCTCCAGCAACACC  
GCAGACCTCTCATCCAGTGTGCTCAAGGAGTTCAAGACCATTCAAGAGGAGGAGTACACAGAGGAGCT  
AGTCACCCAAGGCTGCCCTGATGTTTGAGATCTTAAAAGCAAGCAAGAATGAAGTGATCAGCCAGCAG  
CTGTGCGTCATCTTACACACTGTTACGGTCCCTACCCATCCCCAAGCTCACGGAGATCAAGAGGAAGC  
AGACCTCACGCCTGGACCCTCATTTTCTTAAACAATAAAGAGATGTCCGACGTGACCTTTCTGGTAGAGG  
AAGACCTTCTACGCTCACAAAGTGTCTTGTTCACAGCCTCTCCAGGTTCAAGGCCCTCCTCTCCAGC  
AAGCCAACAAATGACAACACCTGCATAGAGATTGGCTACGTGAAGTATCCCATCTTCCAGCTCGTCATGC  
AGTATCTACTACGGTGGCCAGAGTCACTTCTATTAATAAACAATGAGATCATGGAGCTTCTGTCTGC  
TGCTAAATTTTCCAGCTGGAGGCTCTGCAGCGACACTGTGAGATTATCTGTGCAAAAAGCATCAATACT  
GACAACCTGTGGATATCTACAGCCACGCCAAGTTCCTTGGGGTGACAGAGCTCTCGGCGTACTGTGAAG  
GCTACTTTCTTAAAAATATGATGGTCCTTATTGAAAACGAAGCGTTCAAGCAGCTCCTGTACGACAAAA  
TGGCGAAGGGGCGGCCAGGATGTGCTTACAGACCTACAGAGGACGTTGGCCATTAGGATTCACTATC  
CACCTGTCGTCTTCCAAGGGTCCGTGGTA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR215910 representing NM\_001017525  
 Red=Cloning site Green=Tags(s)


MRKPRPRESREPDGSPVGSGLQRLSHSGVSRPFLVLPPLMEWIRVAVAHAGHRRSFSDSDDDRQAAR  
 LLLPGVDCEPRQLKADDCFCASRKLDAVAIEAKFKQDLGFRMLNCGRTDLVKQAVSLLGPDGINTMSEQG  
 MTPLMYACVRGDEAMVQMLLDAGADLNVEVVSTPHKYPSVHPETRHWALTFAVLHGHIPVVQLLLDAGA  
 KVEGSVEHGEENYSETPLQLAAAVGNFELVSLLLERADPLIGTMYRNGISTTPQGMNSFSQAAAHGHR  
 NVFRKLLAQPEKEKSDILSLEEILAEGTDLAETAPPPLCASRNSKAKLRALREAMYHSAEHGYVDVTIDI  
 RSIGVPWTLHTWLESRLIAFQQHRRPLIQCLLKEFKTIQEEEYTEELVTQGLPLMFEILKASKNEVISQQ  
 LCVIFTHCYGYPYIPKLTEIKRKQTSRLDPHFLNNKEMSDVTFLEGRPFYAHKVLFTASPRFKALLSS  
 KPTNDNTCIEIGYVKYPIFQLVMQYLYGGPESLLIKNNEIMELLSAKFFQLEALQRHCEIICAKSINT  
 DNCVDIYSHAKFLGVTELSAYCEGYFLKNMMVLIENEFKQLLYDKNGEGAGQDVLQDLQRTLAIRIQSI  
 HLSSSKGSVV


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**


**ACCN:** NM\_001017525

**ORF Size:** 1920 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001017525.1](#), [NP\\_001017525.1](#)

RefSeq Size: 3955 bp  
RefSeq ORF: 1923 bp  
Locus ID: 74007  
UniProt ID: [Q6GQW0](#)  
Cytogenetics: 10 C1  
MW: 72.1 kDa