

Product datasheet for **MR215922L4V**

Lace1 (Afg1l) (NM_145743) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lace1 (Afg1l) (NM_145743) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Afg1l
Synonyms:	CG8520; Lace1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_145743
ORF Size:	1443 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215922).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145743.2 , NP_665686.2
RefSeq Size:	2530 bp
RefSeq ORF:	1443 bp
Locus ID:	215951
UniProt ID:	Q3V384
Cytogenetics:	10 22.88 cM



[View online »](#)

Gene Summary:

Putative mitochondrial ATPase. Plays a role in mitochondrial morphology and mitochondrial protein metabolism. Promotes degradation of excess nuclear-encoded complex IV subunits (COX4I1, COX5A and COX6A1) and normal activity of complexes III and IV of the respiratory chain. Mediates mitochondrial translocation of TP53 and its transcription-independent apoptosis in response to genotoxic stress.[UniProtKB/Swiss-Prot Function]