

Product datasheet for **MR215929L3V**

Romo1 (NM_025946) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Romo1 (NM_025946) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Romo1
Synonyms:	2010100O12Rik; AI853864
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_025946
ORF Size:	237 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215929).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025946.6 , NP_080222.1
RefSeq Size:	567 bp
RefSeq ORF:	240 bp
Locus ID:	67067
UniProt ID:	P60603
Cytogenetics:	2 H1
Gene Summary:	Has antibacterial activity against a variety of bacteria including S.aureus, P.aeruginosa and M.tuberculosis. Acts by inducing bacterial membrane breakage (By similarity). [UniProtKB/Swiss-Prot Function]



[View online »](#)