

Product datasheet for MR215943L3

Dusp7 (NM_153459) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Dusp7 (NM_153459) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Dusp7
Synonyms:	AU015694; Mkpx; Pyst2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215943).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_153459
ORF Size:	1266 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_153459.4 , NP_703189.3
RefSeq Size:	3227 bp
RefSeq ORF:	1269 bp
Locus ID:	235584
UniProt ID:	Q91Z46
Cytogenetics:	9 F1
Gene Summary:	Dual specificity protein phosphatase (PubMed:27783954). Shows high activity towards MAPK1/ERK2 (By similarity). Also has lower activity towards MAPK14 and MAPK8 (By similarity). In arrested oocytes, plays a role in meiotic resumption (PubMed:27783954). Promotes nuclear envelope breakdown and activation of the CDK1/Cyclin-B complex in oocytes, probably by dephosphorylating and inactivating the conventional protein kinase C (cPKC) isozyme PRKCB (PubMed:27783954). May also inactivate PRKCA and/or PRKCG (PubMed:27783954). Also important in oocytes for normal chromosome alignment on the metaphase plate and progression to anaphase, where it might regulate activity of the spindle-assembly checkpoint (SAC) complex (PubMed:27783954).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR215943L3