

## Product datasheet for **MR215949L4V**

### **Mocs1 (NM\_020042) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Mocs1 (NM_020042) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mocs1
Synonyms:	3110045D15Rik; AI043170; AW536397
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_020042
ORF Size:	1908 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215949).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_020042.2</a> , <a href="#">NP_064426.2</a>
RefSeq Size:	2640 bp
RefSeq ORF:	1911 bp
Locus ID:	56738
UniProt ID:	<a href="#">Q5RKZ7</a>
Cytogenetics:	17 C



[View online »](#)

**Gene Summary:**

Isoform Mocs1a and isoform Mocs1b probably form a complex that catalyzes the conversion of 5'-GTP to cyclic pyranopterin monophosphate (cPMP). Mocs1a catalyzes the cyclization of GTP to (8S)-3',8-cyclo-7,8-dihydroguanosine 5'-triphosphate and Mocs1b catalyzes the subsequent conversion of (8S)-3',8-cyclo-7,8-dihydroguanosine 5'-triphosphate to cPMP. [UniProtKB/Swiss-Prot Function]