

Product datasheet for **MR215954L3V**

Kctd1 (NM_134112) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Kctd1 (NM_134112) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Kctd1
Synonyms:	4933402K10Rik; AI661543; AW553000
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_134112
ORF Size:	795 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR215954).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_134112.5 , NP_598873.2
RefSeq Size:	1772 bp
RefSeq ORF:	798 bp
Locus ID:	106931
UniProt ID:	Q5M956
Cytogenetics:	18 A1
Gene Summary:	May repress the transcriptional activity of AP-2 family members, including TFAP2A, TFAP2B and TFAP2C to various extent.[UniProtKB/Swiss-Prot Function]


[View online »](#)