

## Product datasheet for **MR216061L3V**

### Azi2 (NM\_001048146) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Azi2 (NM_001048146) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Azi2
Synonyms:	AA410145; AZ2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001048146
ORF Size:	849 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216061).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001048146.1</a>
RefSeq Size:	2579 bp
RefSeq ORF:	852 bp
Locus ID:	27215
UniProt ID:	<a href="#">Q9QYP6</a>
Cytogenetics:	9 F3



[View online »](#)

**Gene Summary:**

Adapter protein which binds TBK1 and IKBKE playing a role in antiviral innate immunity. Activates serine/threonine-protein kinase TBK1 and facilitates its oligomerization. Enhances the phosphorylation of NF-kappa-B p65 subunit RELA by TBK1. Promotes TBK1-induced as well as TNF-alpha or PMA-induced activation of NF-kappa-B. Participates in IFNB promoter activation via TICAM1.[UniProtKB/Swiss-Prot Function]